


miniOne®

S Y S T E M S



CLICK TO SEE SECTION: ELECTROPHORESIS PCR MINILABS REAGENTS INDEX

The Legacy of Winston Walker

Winston Walker dedicated his life to serving his country and finding solutions, whether it be on his aircraft or in industry after his extensive military service. His drive to make science education more accessible led to the development of the MiniOne® Electrophoresis and PCR Systems. His final contribution is The Winston viewer, which helps students of all ages answer various questions around molecules. His impact on science education is global, bringing equality to teachers and students worldwide.

MiniOne® The Winston™

M1050

\$49

See the GLOW! Use fluorescence to detect and see your samples glow in **colors** you won't believe!



Molecules are hard to see with the naked eye—they can be colorless, too small, or not abundant enough. Fluorescence to the rescue! Fluorescence helps you study the invisible. The Winston allows you and your students to investigate phenomena such as: "Does your sample have DNA? Prove it!", "Where is the DNA?", or "What color is chlorophyll?"

Place up to 4 samples on the base and cover with the MiniOne® Photo Hood to reveal the GLOW! Each unit includes:

- One MiniOne® The Winston™ Platform
- One built-in rechargeable battery
- One MiniOne® Photo Hood
- USB charging cable

As a part of the community that Winston helped bring together, we invite you to submit what you discover with The Winston and we will post those protocols online to share with others. For these community-submitted activities check out our website:

<https://www.theminione.com/the-winston>

To submit an activity please email us at: info@minione.com

The Winston will make you **GLOW!**

MiniOne® The Winston™ Platform

M1051

\$34

Includes the base component and charging cable. Great add-on if you already have a MiniOne Photo Hood.



Table of Contents

Electrophoresis Systems.....	4
T-Rack System.....	10
Micropipettes.....	12
Microcentrifuges.....	16
PrepOne™ Sapphire and Photo Hood	17
PCR System.....	18
PCR and Electrophoresis Packages	24
MiniLabs, Electrophoresis.....	28
MiniLabs, PCR.....	36
MiniLab Course Alignment.....	42
GreenGel™ Cups	44
DNA Size Markers	46
Consumables and Plastics.....	48
Resources.....	50
Index	51

Schedule your order to ship when you need it!

Simply let us know future ship dates when you place your order.

Orders must be made out to Embi Tec, and can be placed by:

Email orders@theminione.com

Telephone 1 (858) 684-3190 or 1 (800) 255-1777 (U.S. and Canada)

Fax 1 (858) 684-3195

Online theminione.com (US only)

All prices are in US dollars and are subject to change. We accept PO and credit card payments in US dollars for purchases and shipments to USA addresses. For purchases and shipments to outside of USA, please contact us by phone or email.

MiniOne® Electrophoresis System

M1000-US

\$279


Integrated, real-time electrophoresis system for running DNA labs in one classroom period.



Includes one each of the following items:

(See pages 6 to 9 for complete product details)

- MiniOne® Carriage with blue LED illumination and magnet-activated safety switch
- Gel tank with graphite electrodes
- 42V power supply, 100–240V input
- Amber photo hood for real-time viewing and capturing gel image
- Casting system with casting stand, two gel trays, two reversible combs for six and nine wells, and lid
- **FREE** 2–20 μ L variable volume micropipette (\$59 value!)
- Validation kit: Two GreenGel™ Cups, three DNA samples and TBE buffer concentrate (one per order)



Example results from
Electrophoresis System validation kit
~20 min. run time

MiniOne® Electrophoresis Classroom Package of 10 Systems

M1010-US

~~\$2,790~~ **\$2,499**

One classroom package is for a class of ten student groups, two to three students per group.

Includes:

- Ten sets of MiniOne® Electrophoresis System (**M1000-US**)
- Ten **FREE** 2–20 μ L variable volume micropipette (\$590 value!)
See new Micropipette list price
- Validation kit: Two GreenGel™ Cups, three DNA samples and TBE buffer concentrate (one per order)




MiniOne® Carriage

M2007

\$249

Housing and control unit with sealed LED lights and power controls.

- + and – electrical contacts for carbon electrodes of gel tank
- Design ensures correct orientation fit
- Two rows of blue LEDs illuminate the gel from the sides
- Choose from two intensities—low light for loading, bright light for DNA viewing
- Built-in, magnet-activated on/off power switch




MiniOne® Photo Hood

M2005

\$29

Molded photo hood for 360° viewing with ventilation slots and black imaging platform.

- Four embedded magnets activate power carriage
- Amber filter blocks blue light and transmits green light showing bright fluorescent DNA bands against a dark background
- Ventilation slots prevent condensation inside the hood



NEW!

MiniOne® Photo Hood Phone Platform

M2017 (Set of 5)

\$25

Molded adapter allows you to place a phone on top of the MiniOne Photo Hood to easily take videos, time lapses, or still images of your electrophoresis run. (Photo Hood not included)

- Latches into the photo hood vents
- Can support mobile devices up to 6.2" (L) x 3.0" (W)" long
- Hands-free documentation



MiniOne® 42V Power Supply

M2006

\$19

- Input voltage: 100–240 V, 50/60 Hz, 0.5A
- Output voltage: 42V, 0.19A
- Certifications: CE, FCC, cULus, RoHS compliant




MiniOne® Gel Tank

M2001

\$49

Molded polycarbonate tank with graphite electrodes.

- Three ridges on walls of tank to guide correct orientation of gel tray
- Durable 1 cm graphite electrodes replace standard fragile platinum wires for easy cleaning
- Black and silver gray gel platform provides contrasting background for imaging gels




MiniOne® Gel Tray Platforms

M2014 Black gel tray platform (10 per pack)


M2015 Silver gray gel tray platform (10 per pack)

\$19 per pack

Slots on the sides of the platform ensure that it only fits into the tank in the correct orientation.



Use with DNA separation. Black color plastic plate printed with a fluorescent green ruler on the edge, texture on well area.



Use with color dyes separation. Silver gray color plastic plate with texture on well area.


MiniOne® Gel Casting System

M2002

\$49

A compact and self-contained unit with two clear gel trays and two gel combs that fit inside a casting stand with lid.

- Molded casting stand eliminates leakage—no tape needed
- Gel tray guides create a one-way fit for the gel trays when casting, ensuring the proper positioning of the wells for use in the gel tank
- Dual reversible comb with one beveled side, 1.5 mm max. thickness, six and nine wells, and white lid for overnight storage of pre-poured gels
- Beveled comb for easy comb removal and sample loading
- White cover blocks light so you can keep the gels you pour ahead of time safe from light exposure



MiniOne® Gel Trays

M2013 (10 per pack)

\$49

M2013-Bulk (50 per pack)

\$195

New gel tray design, clear acrylic trays with one orientation fit feature.

- Use in MiniOne tank and casting stand



NEW!

MiniOne® Casting Stand Cover

M2018 (5 per pack)

\$34

- Fits MiniOne Gel Casting Stand
- Blocks light for advanced gel prep




MiniOne® Gel Combs

M2004 (Set of 2)

\$19

- Dual reversible gel comb with one beveled side
- 1.5 mm max. thickness, 6+6 / 9+9 wells



NEW!

Colorful MiniOne® T-Rack™ System

The modular MiniOne T-Racks are designed to give you the most flexibility with how you set up your lab stations.

MiniOne® T-Rack™ MicrotubeRack - 1.5 and 2.0 mL

M3181 (Set of 5 racks)

\$49

Five assorted colors, holds 3 x 4 of 1.5 or 2.0 mL microtubes per rack. (Tubes not included)



Capacity and configuration guide



Colors may vary

MiniOne® T-Rack™ PCR Tube Rack - 0.2 and 0.5 mL

M3180 (Set of 5 racks)

\$49

Five assorted colors, each rack holds 2 x 5 of 0.2 mL PCR tubes and 3 x 4 of 0.5 mL microtubes. (Tubes not included)



Capacity and configuration guide



Colors may vary

MiniOne® T-Rack™ Micropipette Tip Rack and Cover

M3139 (Set of 5 racks)

\$49

Five assorted rack colors with clear lids, pre-racked with 2–200 µL Universal tips, low binding, 36 tips/rack.



Capacity and configuration guide



Colors may vary

MiniOne® T-Rack™ System Tray

M3182 (Set of 5 trays)

\$15

Trays for holding up to three T-Rack modules per platform for easy workstation organization.




MiniOne® T-Rack Combo Pack

M3143 (Set of 5)

~~\$163~~ **\$149**

Assorted colors (5), complete set of racks for 0.2 mL/0.5 mL PCR tubes, 1.5/2.0 mL microtubes, and 2–200 µL universal tips, with convenient tray to hold up to 3 racks at a time. Combo pack contains one set each of M3139, M3180, M3181, M3182, to make 5 complete sets.



NEW!

MiniOne® Mini Erlenmeyer Flasks

M2019 (Set of 5)

\$19

Safer than traditional glass, this affordable set of 5 polypropylene 170 mL Erlenmeyer flasks is ideal for the biology classroom. They're light, rugged, autoclavable, and feature a pour spout, making them perfect for storing and pouring buffer. Minimize spills and breakage in the lab with MiniOne Mini Erlenmeyer Flasks. Non-microwavable, do not freeze.

- 170 mL (6 oz) graduated every 25 mL (1 oz)
- Polypropylene
- Autoclavable (*remove sticker before autoclaving*)



MiniOne® Micropipettes

M2008, M2010, M2011, M2012

Best Classroom Pipette at the Best Price!

\$59

One micropipette per pack; comes with inspection certificate, calibration tool, operation manual and sample pipette tips.



Plunger button

- Low pipetting force
- Color-coded by volume range

Volume adjuster

- Smooth volume adjustment on plunger
- Easy-turn dial

Ejector button

- Low operating force

Volume display

- 3 digit display
- Decimal place is color-coded in red

Universal tip holder

- Accepts most leading brand pipette tips

Features:

- Adjustable volume micropipette with durable and reliable quality
- Ergonomic design provides comfortable operation for small or large hands
- Two 'stops' on the plunger to allow for accurate and complete sample dispensation
- Ultra-affordable for educational use
- Can be self-calibrated

MiniOne® Micropipette Set

M2016

\$195

This set of three lab-quality micropipettes comes protected in a case constructed of durable, easy-to-clean, waterproof-coated Lycra, ensuring your pipette set stays safe, organized, and clean.



Includes one of each variable volume micropipettes:

M2008 2–20 µL

M2010 20–200 µL

M2011 100–1000 µL

MiniOne® Micropipette Stand

M2021

\$75

Keep your bench organized!

Holds up to nine MiniOne Micropipettes.



Cat. No.	MiniOne® Micropipette Volume Range	Accuracy	Repeatability	Price
M2008	2–20 µL (H20)	± 2.5-1.0%	≤ 1.50–0.30%	\$59.00
M2010	20–200 µL (H200)	± 1.8-0.6%	≤ 0.50–0.15%	\$59.00
M2011	100–1000 µL (H1000)	± 1.5-0.6%	≤ 0.30–0.15%	\$59.00
M2012	1–10 µL (H10)	± 2.5-1.0%	≤ 1.50–0.40%	\$59.00
M2016	Set of 3 MiniOne® Micropipettes, 2–20 µL, 20–200 µL and 100–1000 µL adjustable volume micropipettes (one each) in a protective carrying case			\$195.00
M2021	Pipette stand for the MiniOne® micropipettes; hold 9 micropipettes			\$75.00

All micropipettes are manufactured according to the standards ISO13485. QC inspection and test are complying with standard ISO8655.


Micropipette Tips - See page 49 for specifications		
M3112	1-10 µL, pk of 250 tips 1-10 µL, pk of 250 tips	\$9.50
M3111	2–200 µL, pk of 250 tips	\$9.50
M3134	2–200 µL, pk of 1,000 tips	\$28.50
M3118	100–1,000 µL, pk of 250 tips	\$9.50
M3139	Micropipette Tip Rack and Cover, set of 5 assorted colors • 36 x 2–200 µL Universal tip per rack • See page 10 for details	\$49.00

ONE Series™ Micropipettes

EA-1001 to EA-1006

\$89

One micropipette per pack; comes with inspection certificate, calibration tool, operation manual and sample pipette tips.



Plunger button

- Low pipetting force
- Color-coded by volume range

High resting hook

- Short thumb reach
- Comfortable for small and large hands

Balanced weight

- Minimizes fatigue
- Optimal weight distribution

Volume adjuster

- Smooth, tactile clicks
- Easy-turn dial

Ejector button

- Low operating force

Volume display

- 4 digit display
- Color-coded by volume range

Universal tip holder

- Accepts most leading brand pipette tips

Features:

- Adjustable volume micropipette combines all the features of accuracy, reliability, and ergonomic design in one
- Great value for scientific research and educational uses
- Autoclavable
- Can be self-calibrated

ONE Series™ Micropipette Stand

EA-1011

\$75

Keep your bench organized!
Holds up to nine ONE Series Micropipettes.



Cat. No.	ONE Series™ Micropipette Volume Range	Accuracy	Repeatability	Price
EA-1001	0.1–2 µL	± 12.0–1.5%	≤ 6.00–0.70%	\$89.00
EA-1002	2–20 µL	± 2.5–1.0%	≤ 1.50–0.30%	\$89.00
EA-1003	20–200 µL	± 1.8–0.6%	≤ 0.50–0.15%	\$89.00
EA-1004	10–100 µL	± 1.8–0.8%	≤ 0.50–0.15%	\$89.00
EA-1005	100–1000 µL	± 1.5–0.6%	≤ 0.30–0.15%	\$89.00
EA-1006	0.5–10 µL	± 2.5–1.0%	≤ 1.50–0.40%	\$89.00
EA-1010	Set of 3 ONE Series Micropipettes: 2–20 µL, 20–200 µL and 100–1000 µL adjustable volume micropipettes and MiniOne pipet holder			\$265.00
EA-1011	Pipette Stand for the ONE Series Micropipettes; holds 9 micropipettes			\$75.00
EA-1014	Set of 4 ONE Series Micropipettes: 0.1-2 µL, 2–20 µL, 20–200 µL and 100–1000 µL adjustable volume micropipettes (one each), w/Holder			\$349.00

All micropipettes are manufactured according to the standards ISO13485. QC inspection and test are complying with standard ISO8655.

Micropipette Tips - See page 49 for specifications		
M3112	1-10 µL, pk of 250 tips 1-10 µL, pk of 250 tips	\$9.50
M3111	2–200 µL, pk of 250 tips	\$9.50
M3134	2–200 µL, pk of 1,000 tips	\$28.50
M3118	100–1,000 µL, pk of 250 tips	\$9.50
M3136	2–200 µL, ten racks • Tips packed in racks of 96 micropipette tips for easy use and storage • Has all features of the bulk packaging version	\$55.00

MiniOne® Centrifuges

Portable mini centrifuges ideal for bringing small droplets to the bottom of tubes, for micro-filtrations, or basic separations. Brushless motor for low noise level. CE marked.

MiniOne® Multi Speed Centrifuge

M2031

\$350

Features:

- Adjustable speed from 1,000 to 10,000 RPM
- One combi-rotor with positions for varied sizes of microcentrifuge tubes and PCR tubes—no need to change rotors



Specifications	
Input Voltage	100-240V AC, 50-60Hz
Rotor Capacity	2 x 8 x 0.2 mL PCR tubes and 6 x 1.5/2.0 mL and 6 x 0.5 mL microcentrifuge tubes
Speed Range	Adjustable speed: 1,000 to 10,000 RPM, max 3,200 x g
Timer	15 sec to 99 minutes or continuous operation
Dimensions	20 x 16 x 13 cm
Weight	1.05 kg

MiniOne® Single Speed Centrifuge

M2032

\$150

Features:

- Easy operation: one button to control ON or OFF
- Fixed speed at 10K RPM providing 4,800 x g (RCF) when using the microcentrifuge tube rotor



Specifications	
Input Voltage	100-240V AC, 50-60Hz
Rotor Capacity	PCR tube rotor: 2 x 8 x 0.2 mL PCR tubes, Microcentrifuge tube rotor: 6 x 1.5/2.0 mL or 6 x 0.5 mL with adaptors or 6 x 0.2 mL with adaptors
Speed Range	Fixed speed, 10,000 RPM
Dimensions	17 x 15 x 13 cm
Weight	0.95 kg

PrepOne™ Sapphire and Photo Hood

M5000

\$549

A blue LED illuminator for smart devices to take gel images with ease. Adapt your current electrophoresis system for safe blue light illumination and non-toxic GelGreen™ stain.

Package includes:

- One PrepOne™ Sapphire Blue LED Illuminator
- One photo hood for smart devices (assembly required)
- One amber filter

Photo hood


- Compact and collapsible design turns your bench into an instant dark room.
- Wide top platform holds most smart devices
- Height: 15.4 cm, provides a suitable focal length for most phone cameras
- Base DIM: 13.8 (W) x 13.8 (L) cm

Amber filter

- Cuts out excess blue light to enhance the contrast of the DNA signal
- DIM: 15.2 (W) x 15.2 (L) x 0.3 (H) cm

PrepOne™ Sapphire

- Side illuminating blue light allows direct visualization of DNA bands
- 4x4 blue LEDs last over 50,000 hours
- Two built-in light intensities
- DIM: 21.4 (W) x 16.3 (L) x 1.5 (H) cm





MiniOne® PCR System


M4000-US

Teach and Do PCR Labs in 45 minutes!


\$799

Fast PCR thermal cycler controlled by App on Bluetooth® LE mobile device (see page 22). Amplify DNA in one classroom period.

- One MiniOne thermal cycler for PCR with 16 sample wells
- One 12V, 100W power supply, 100–240 VAC
- One **FREE** validation kit
- Temperature range **4°C (below ambient) to 99°C**
- Ability to pause run to evaluate cycle number



CE marked



Lane 1: MiniOne DNA marker (see page 34)
 Lane 2: PCR fragment 1
 Lane 3: PCR fragment 2
 Lane 4: PCR fragment 3
 Lane 5: negative control

Example results from PCR System validation kit ~ 20 min. amplification and ~ 20 min. run time

Validation kit provides enough reagents for three PCR runs, each amplifying three different size fragments, then visualizing the results with MiniOne® Electrophoresis

Features:

- Peltier cooling technology and custom algorithm drives **fast thermal** cycling
- **Mobile App** for programming and monitoring via **Bluetooth® LE**
- Indicator lights on front show status of the run
- Fully compatible with standard reagents, consumables, and protocols

Benefits:

- Complete a PCR protocol in a single class period
- Intuitive programming interface for students
- Keep your samples at 4°C when done



Specifications

Sample capacity	16 x 0.2 mL standard PCR tubes
Temperature range	4°–99°C
Heated lid	Yes, with safety interlock
Communications	Bluetooth® Low Energy wireless technology
Software	Graphical programming interface; Real-time protocol monitoring
Weight	1.9 lb (860 g) approx.
Dimensions	12 x 12 x 12 cm (4.7 x 4.7 x 4.7 in) approx.
Operating voltage	100–240 VAC



MiniOne® PCR System (continued)

16 sample capacity

Perfect for two to four student groups or a small class.

For larger classes, add more systems to maintain student exposure and participation.



Active heating and cooling

The Peltier element drives rapid temperature transitions and faster protocols. Keeps your samples at 4°C at the end of the run.

Specialized algorithm

Custom control algorithm is the key for fast cycling and precise thermal control.

Visual indicators

Bright LEDs show the state of the machine and progress of the protocol.



Compact

The MiniOne® PCR System takes up minimal space on your bench, stores easily, and can be moved wherever it's needed.



Safe

Fully enclosed system for safe operation in the classroom.




Heated Lid

No condensation or messy mineral oil. Safety switch keeps protocol from starting until the lid is closed.

How fast is it?

A typical fast protocol - MiniOne® PCR System is done in **50%** of the time.

Step	Temperature	Duration	Cycles
Denaturation	94°C	5 sec	20 cycles
Annealing	54°C	5 sec	
Extension	72°C	5 sec	




MiniOne® PCR App

Intuitive, student-centered app for programming and monitoring your PCR protocols.

Features:


- Interactive screens lead students through each step of setting up the PCR protocol
- **Pause** feature stops the protocol at the end of the extension step for convenient cycle number analysis
- Graphical output **displays real-time temperature data** that can be saved or emailed at the end of the run
- **Linked protocol** to automatically run consecutive protocols



iPhone (iOS 9 or higher)



Android phone (Lollipop or higher)



iPad/Mini iPad



MiniOne® PCR



Android controller



MacBook or PC Laptop*



Chromebook (Bluetooth 4.1 or above)

*Requires MiniOne® Bluetooth® LE Dongle **M4060 \$99**



Android Mobile Controller

M4050

\$99

Android mobile controller with MiniOne® PCR App pre-loaded.

- 7-inch, full color graphic display for a modern, student-centered interface
- Bluetooth® LE connectivity for programming and monitoring your MiniOne® PCR System



1
Bluetooth® LE connection

View and connect to available PCR Systems with a wireless Bluetooth® LE connection. No cables or wifi needed!

2
Setup

Create a new protocol, access your last protocol or browse the library of saved protocols from the Setup menu.

3
Protocol

Select constant temperature mode for restriction digests, transformation and DNA extractions, or PCR mode for DNA amplification.

4
Constant Temperature

Set up constant temperature protocols to incubate samples at 0-99°C using a simple, intuitive interface.

5
PCR

Graphic PCR programming interface allows students to enter time and temperatures directly on the PCR overview graph.

6
Cycling

The run screen shows a continuously updated, real-time readout of the temperature, current cycle, current step, and elapsed time.

MiniOne® PCR with Android Mobile Controller

M4001

\$898

- One MiniOne® PCR System with one validation kit - **CE** marked
- One Android Mobile Controller with MiniOne® PCR App installed



MiniOne® PCR/Electrophoresis Package I

M4011

\$995

- One MiniOne® PCR System with one validation kit - **CE** marked
- One MiniOne® Electrophoresis System - **CE** marked
- One **FREE** 2–20 μ L variable volume micropipette

For a group of two to three students



Package II

M4012

\$1,259

- One MiniOne® PCR System with one validation kit - **CE** marked
- Two MiniOne® Electrophoresis Systems - **CE** marked
- Two **FREE** 2–20 μ L variable volume micropipettes

For two groups of students



Package III

M4026

\$2,995

- Two MiniOne® PCR Systems with two validation kits - **CE** marked
- Six MiniOne® Electrophoresis Systems - **CE** marked
- Six **FREE** 2–20 μ L variable volume micropipettes

For six groups of students





MiniOne® PCR/Electrophoresis Package IV

M4039

\$4,475

- Three MiniOne® PCR Systems - **CE** marked
- Two PCR validation kits
- Nine MiniOne® Electrophoresis Systems - **CE** marked
- Nine **FREE** 2–20 µL variable volume micropipettes

For nine groups of students



MiniOne® Bioscience Classroom Starter Package V

M4040

\$4,990

This bioscience classroom starter package includes the essential equipment for teaching hands-on molecular biology concepts using DNA amplification and separation in the classroom. The starter package is ideal for nine groups of two to three students, and includes:

- Three MiniOne® PCR Systems - **CE** marked
- Two PCR validation kits
- Nine MiniOne® Electrophoresis Systems - **CE** marked
- Three MiniOne® Single Speed Microcentrifuges - **CE** marked
- One MiniOne® Micropipette Set with Case, includes one of each: 2–20, 20–200, 100–1000 µL variable volume micropipettes
- Nine **FREE** 2–20 µL variable volume micropipettes

For nine groups of students

Includes micropipette set for teacher prep!



MiniOne® MiniLabs

Our hands-on MiniLabs are a fun and engaging series of modules that take students from mastery of basic biotech skills, through popular applications of electrophoresis in forensics, DNA fingerprinting, and human genetics, and finally, to a challenging, real-world investigation of a foodborne outbreak.

Gel Loading Practice MiniLab

M3002

\$49



Teach your students how to read, adjust, and use a micropipette. Practice pipetting and loading samples into the wells of a real gel before handling valuable DNA samples.

For middle school students and any first-time micropipette users (grades 4–12).

Each MiniLab contains enough materials for 10 workstations, 2–3 students per workstation.


Materials include:

Twenty precast agar plates with 3 rows of 8 wells per row	One bag of 1.7 mL microcentrifuge tubes
Twenty adhesive practice pipetting sheets	One bag of 2–200 µL micropipette tips
Two tubes of dye samples (10 mL each)	Teacher's Guide

Laminated Practice Pipette Cards

M2022 (Set of 20)

\$24



Reusable practice pipette cards from our Gel Loading Practice MiniLab. Lamination makes it easy to clean so you are ready for your back-to-back classes.

Colorful Dye Electrophoresis MiniLab

M3007

\$59



This fun and colorful lab introduces students to gel electrophoresis principles, including the basics of electricity and macromolecules. Predict how molecular size and electrical charge affect a molecule's migration in a separation matrix.

Appropriate for middle school and beginning high school students (grades 7–10)

Each MiniLab contains enough materials for 10 workstations, 2–3 students per workstation.

Materials include:

Ten 1% agarose gel cups	One bag of 0.65 mL microcentrifuge tubes
Nine color dye samples	One bag of 2–200 µL micropipette tips
One bottle of 100 mL Tris-Borate-EDTA (TBE) buffer concentrate	Teacher's Guide

Each MiniLab is an all-inclusive science lab kit. The ready-to-pour gel cups greatly simplify steps to make agarose gels. The TBE buffer concentrate and easy-to-load DNA/color dye samples assure good results, minimize teacher prep work, and maximize student success in the lab.

Candy Color Electrophoresis MiniLab

M3009

\$59



Investigate the phenomenon of food dyes using candy to illustrate the effect of mass and charge during electrophoresis.

Appropriate for middle school and beginning high school students (grades 7–10).

Each MiniLab contains enough materials for 10 workstations, 2–3 students per workstation.

Materials include:

Ten 1% agarose gel cups	Dye extraction buffer
One bottle of 100 mL Tris-Borate-EDTA (TBE) buffer concentrate	Candies in six colors
One bag of 1.7 mL microcentrifuge tubes	Ten 10-well dye extraction trays
One bag of 2–200 µL micropipette tips	Teacher's Guide

Laminated Gel Annotation and Photo Template

M2023 (Set of 10)

\$14



Place your gel on these laminated cards make the color dyes POP after your electrophoresis run! Dual sided for 6- or 9-well dye electrophoresis activities.

Determining the Genetics of a Ca\$H Cow MiniLab

M3011

\$59



Using an understanding behind cheese production, genetic inheritance, and gel electrophoresis, students will determine genotype of two bulls and three cows, and recommend which combination a dairy farmer should purchase to produce more valuable offspring.

Appropriate for middle school life science classes, genetics classes, and agriculture courses (grades 7–12).

Each MiniLab contains enough materials for 10 workstations, 2–3 students per workstation.

Materials include:

Eight color dye samples	One bag of 0.65 mL microcentrifuge tubes
Ten 1.5 % agarose gel cups	One bag of 2–200 µL micropipette tips
One bottle of 100 mL Tris-Borate-EDTA (TBE) buffer concentrate	Teacher's Guide

MiniOne® MiniLabs (continued)

Electrophoresis 101 MiniLab

M3001

\$39


Students are introduced to the principles of gel electrophoresis by separating colorful dyes and DNA samples on an agarose gel. Challenge your students' analytical and mathematical skills as they construct a standard curve to determine the sizes of unknown DNA fragments.

Appropriate for high school students (grades 9–12).

Each MiniLab contains enough materials for 10 workstations, 2–3 students per workstation.

Materials include:

Ten 1% agarose GreenGel™ Cups
 Three color dye samples and three DNA samples
 One bottle of 100 mL Tris-Borate-EDTA (TBE) buffer concentrate
 One bag of 0.65 mL microcentrifuge tubes
 One bag of 2–200 µL micropipette tips
 Teacher's Guide



DNA Fingerprinting MiniLab

M3004

\$69


How is DNA used to trace the history and heritage of an individual? Students help scientists identify the father of a baby humpback whale using DNA fingerprinting technology. Engage your students with a real-world application of genetics as they analyze a complex array of DNA bands to arrive at a logical solution.

Appropriate for high school students (grades 9–12).

Each MiniLab contains enough materials for 10 workstations, 2–3 students per workstation.

Materials include:

Ten 1% agarose GreenGel™ Cups
 Five DNA samples
 One bottle of 100 mL Tris-Borate-EDTA (TBE) buffer concentrate
 One bag of 0.65 mL microcentrifuge tubes
 One bag of 2–200 µL micropipette tips
 Teacher's Guide



CSI Forensics MiniLab

M3005

\$69


Explore a crime scene investigation in the classroom. Learn an exciting real-world application of gel electrophoresis and the statistical principles of human genetic identification. Students will logically integrate multiple lines of evidence, including fingerprints, hair samples, and DNA fingerprinting to connect an individual to a crime scene and solve the mystery of "Who Killed Dr. Ward?"

Appropriate for high school students (grades 9–12).

Each MiniLab contains enough materials for 10 workstations, 2–3 students per workstation.

Materials include:

Ten 1% agarose GreenGel™ Cups
 Five DNA samples
 One bottle of 100 mL Tris-Borate-EDTA (TBE) buffer concentrate
 One bag of 0.65 mL microcentrifuge tubes
 One bag of 2–200 µL micropipette tips
 Teacher's Guide



PTC Genetics MiniLab

M3003


\$69

Solve a genetic mystery using gel electrophoresis. Students develop a hypothesis about the inheritance of a trait in a family, then test their hypotheses by running restriction fragments on a gel, analyzing a Punnett square, and constructing a family tree.

Appropriate for high school students (grades 9–12).

Each MiniLab contains enough materials for 10 workstations, 2–3 students per workstation. Materials include:

Ten 2% agarose GreenGel™ Cups
 Six pre-digested DNA samples
 Forty pieces of PTC tasting and taste control papers
 One bottle of 100 mL TBE buffer concentrate
 One bag of 0.65 mL microcentrifuge tubes
 One bag of 2–200 µL micropipette tips
 Teacher's Guide



MiniOne® MiniLabs (continued)

PTC Inheritance and Graphical Analysis MiniLab

M3012

\$99

Explore Mendelian genetic inheritance, use Punnett Squares to make a claim, and see your evidence with DNA electrophoresis. This lab also features DNA fragment size analysis and demystifies why the standard curve is plotted as Log-Y!


Appropriate for AP Biology, Honors and Advanced Biology (grades 9-12).

Each MiniLab contains enough materials for 10 workstations, 2–3 students per workstation.

Materials include:

Ten 2% agarose GreenGel™ Cups
Six pre-digested DNA samples
MiniOne® DNA Marker
Forty pieces of PTC tasting and taste control papers

One bottle of 100 mL Tris-Borate-EDTA (TBE) buffer concentrate
One bag of 0.65 mL microcentrifuge tubes
One bag of 2–200 µL micropipette tips
Teacher's Guide



Hunting the Inheritance of Huntington's Disease MiniLab

M3010

\$99

In this lab, students will examine family history to construct a pedigree and will assess molecular data to make predictions about inheritance of the disease in fraternal twins. Students will then perform gel electrophoresis to confirm genotype.

Appropriate for high school life science, genetics classes and AP Biology.

Each MiniLab contains enough materials for 10 workstations, 2–3 students per workstation.

Materials include:

Ten 2% agarose GreenGel™ Cups
Four Ready-to-Load DNA samples
100 bp DNA Ladder
One bottle of 100 mL Tris-Borate-EDTA (TBE) buffer concentrate
One bag of 0.65 mL microcentrifuge tubes
One bag of 2–200 µL micropipette tips
Teacher's Guide



Restriction Digest Basics MiniLab

M6050

\$69

Cells have mechanisms for cutting long strands of nucleic acid into shorter strands—a type of molecular scissors. There are several reasons why cells need to cut their DNA or RNA. In this lab students will explore what restriction enzymes do and determine electrophoresis fragment sizes by comparing bands to a molecular weight standard.


Appropriate for high school students (grades 9-12).

Each MiniLab contains enough materials for 10 workstations, 2-3 students per workstation.

Materials include:

Ten 1.5% agarose GreenGel™ Cups
Four DNA samples (three pre-digested, one undigested)
MiniOne® Universal DNA Marker

One bottle of 100 mL Tris-Borate-EDTA (TBE) buffer concentrate
One bag of 0.65 mL microcentrifuge tubes
One bag of 2–200 µL micropipette tips
Teacher's Guide



Restriction Analysis of DNA MiniLab

M6053

\$120

“Molecular scissors” for cutting DNA or RNA are used by cells for many reasons, including defense against an invading host, or genetic recombination. These restriction enzymes allow researchers to study smaller pieces of DNA more thoroughly. In this restriction digestion lab, students will explore what restriction enzymes do, perform a single and double digest of a synthetic piece of DNA, predict fragment sizes and compare to fragment size of the actual digested DNA run on an agarose gel.


Appropriate for AP Biology, Honors and Advanced Biology (grades 9-12).

Each MiniLab contains enough materials for 10 workstations, 2–3 students per workstation.

Materials include:

Ten 1.5% agarose GreenGel™ Cups
Undigested DNA Sample
Four pre-digested DNA controls
MiniOne® Universal DNA marker
Enzyme dilution buffer
Two restriction enzymes
MiniOne® 5X Sample Loading Dye

One bottle of 100 mL Tris-Borate-EDTA (TBE) buffer concentrate
One bag of 0.65 mL microcentrifuge tubes
One bag of 2–200 µL micropipette tips
One bag of 0.2 mL PCR tubes
Teacher's Guide



MiniOne® MiniLabs (continued)

Foodborne Outbreak Investigation MiniLab

M3006

\$120


This MiniLab is a student-driven discovery process based on a real *Shigella* outbreak in 2000. Students use scientific reasoning and forensic science principles to analyze epidemiological data, develop a hypothesis, and test their hypothesis with gel electrophoresis. They will tabulate data expressed as text to systematically analyze the case and evaluate experimental approaches used by their team and others in the class. Along the way they will develop an understanding of foodborne outbreaks and foodborne illness, topics that are frequently in the news and relevant to students' lives.

Appropriate for high school students (grades 9–12), AP, honors, and advanced biology students.

Each MiniLab contains enough materials for 10 workstations, 2–3 students per workstation.

Materials include:

- Ten 1% agarose GreenGel™ Cups
- Eleven DNA samples
- One bottle of 100 mL Tris-Borate-EDTA (TBE) buffer concentrate
- Two bags of 0.65 mL microcentrifuge tubes
- One bag of 2–200 µL micropipette tips
- Teacher's Guide
- A two-part detailed guide with background information, step-by-step procedures and worksheets for students
- Quiz question bank and answers



	1	3	7	9	10	11	2	4	5	6	8	12
Potato Chips		X	X	X		X	X	X				
French Onion Dip	X	X	X	X	X	X	X	X		X	X	X
Tortilla Chips		X	X	X								X
Salsa	X							X			X	X
Guacamole				X					X			X
Homemade Potato Salad		X	X	X	X	X	X	X	X	X		X
Burgers						X	X	X				X
Hot Dogs		X	X	X	X	X	X	X				X
Five Layer Bean Dip		X	X	X	X	X	X	X				X

NGSS-Aligned Color Dyes and Gel Electrophoresis MiniLab

M3008

\$199


A comprehensive, 5E inquiry, week-long lesson plan that introduces students to separation science in a variety of contexts. Student-driven inquiry activities impart the scientific background needed to understand gel electrophoresis. The curriculum culminates in a gel electrophoresis experiment using colorful dyes. The curriculum is provided as a downloadable PDF manual which includes comprehensive background readings, student worksheets, and teacher support materials. Lab materials and a PowerPoint presentation are also included.

Recommended for middle school students (grades 6-9).


Each MiniLab contains enough materials for 10 workstations, 2-3 students per workstation.

Materials include:

- Ten 1% agarose gel cups
- Six color dye samples
- One set of filter papers and green food coloring for exploring paper chromatography
- One bottle of 100 mL Tris-Borate-EDTA (TBE) buffer concentrate
- One bag of 0.65 mL microcentrifuge tubes
- One bag of 2–200 µL micropipette tips
- One copy of the curriculum



10 minutes run time



Start

Finish

Bundle and Save!

M3008 3-Pack Reagents for 30 workstations, 1 copy of the curriculum	\$337
M3008 5-Pack – Reagents for 50 workstations, 1 copy of the curriculum	\$475



MiniOne® MiniLabs (continued)

PCR 101 MiniLab: Amplification from the Lambda Phage Genome

M6001

\$99



In this hands-on PCR MiniLab, students use polymerase chain reaction (PCR) to amplify three segments of the Lambda phage genome. They will look at sequence data, predict the fragment sizes of the PCR products, then compare their predictions to the PCR products they amplify and run on an agarose gel. Complete amplification in 17 minutes with the MiniOne® PCR system and the FastTaq™ Master Mix.

Appropriate for high school students (grades 9–12), AP, honors, and advanced biology students.

Each MiniLab contains enough materials for 10 workstations, 2–3 students per workstation. **Materials include:**

- | | |
|--|--|
| Ten 2% agarose GreenGel™ Cups | MiniOne® 5X Sample Loading Dye |
| FastTaq™ PCR MasterMix (2X) | One bag of 0.2 mL thin-wall PCR tubes |
| Three primer sets, both forward and reverse primers included in each set | One bag of 0.65 mL microcentrifuge tubes |
| Lambda phage genomic DNA | One bottle of 100 mL Tris-Borate-EDTA (TBE) buffer concentrate |
| Sterile nuclease-free water | Teacher's Guide |
| MiniOne® DNA marker | |

Bundle and Save!


COMPLETE PACKAGE!

	M6001 \$99	M6002 \$85	M6004 \$129	M6003 \$139
M6002 - PCR 101 MiniLab, PCR Reagents (no electrophoresis reagents)	✓	✓	✓	✓
M3103TBE - Ten 2% agarose GreenGel™ Cups with TBE buffer concentrate (see page 45)	✓			✓
M3136 - 10 racks of micropipette tips (2–200 µL) (see page 49)			✓	✓

PCR Cycle Number Analysis MiniLab

M6005

\$99




Visualize the power of exponential growth with PCR! Students will set up PCR reactions and analyze the products after a variable number of cycles. Students will estimate the minimum number of cycles needed to detect a PCR product on an agarose gel, and visualize and appreciate exponential growth.

Appropriate for high school students (grades 9–12), AP, honors, and advanced biology students.

Each MiniLab contains enough materials for 10 workstations, 2–3 students per workstation. **Materials include:**

- Ten 1% agarose GreenGel™ Cups
- MiniOne® DNA marker
- MiniOne® 5X Sample Loading Dye
- FastTaq™ PCR MasterMix (2X)
- One primer set, forward and reverse
- Lambda phage genomic DNA
- One bag of 0.65 mL microcentrifuge tubes
- One bag of 0.2 mL thin-wall PCR tubes
- One bottle of 100 mL TBE buffer concentrate
- Teacher's Guide



Need some racked tips? or PCR tube racks?



Try the MiniOne® T-Rack Micropipette Tip Rack and Cover (**M3139**) or the MiniOne T-Rack™ PCR Tube Rack - 0.2 and 0.5 mL (**M3180**) (see pg. 10 for details)

MiniOne® MiniLabs (continued)

A Taste of Genetics MiniLab: Extract and Amplify the PTC Gene

M6010

\$115



This hands-on MiniLab introduces students to the science of human genetic variation through DNA extraction, PCR amplification, restriction digest, and analysis of the TAS2R38 taster gene of their own DNA, and compare genotype to phenotype.

Appropriate for high school biology students, especially honors and advanced placement, and college level biology.

Each MiniLab contains enough materials for 10 workstations, 2–3 students per workstation.

Materials include:

- | | |
|---|--|
| DNA extraction solution | One bag of 0.65 mL microcentrifuge tubes |
| Forward and reverse primers for PTC genes | Forty pieces of PTC taste paper and taste control papers |
| Taq polymerase master mix (2X) | Ten 2% agarose GreenGel™ Cups |
| HaeIII restriction enzyme | One bottle of 100 mL Tris-Borate-EDTA (TBE) buffer concentrate |
| Restriction enzyme dilution buffer | Two grams table salt |
| MiniOne® 5X Sample Loading Dye | Teacher's Guide |
| MiniOne® DNA Marker | |
| One bag of 0.2mL thin-walled PCR tubes | |

Bundle and Save!

	M6010 \$115	M6012 \$99	M6013 \$155
M6012 - A Taste of Genetics MiniLab, DNA Extraction and PCR Reagents (no electrophoresis reagents)	✓	✓	✓
M3103TBE - Ten 2% agarose GreenGel™ Cups with TBE buffer concentrate (see page 45)	✓		✓
M3136 - 10 racks of micropipette tips (2–200 µL) (see page 49)			✓

COMPLETE PACKAGE!

Extension Activities

available at <https://theminiOne.com/free-classroom-activities/>

Introduction to NCBI Bioinformatics – Students will learn to navigate NCBI, explore the different types of information available, and apply it to look deeper at the TAS2R38 gene, including sequence, chromosomal location, and associated phenotypes and variants.

Hardy-Weinberg – Students will create a simulation of allele frequencies in a population using principles from the Hardy-Weinberg model and a TI-84 Plus graphing calculator or computer including evaluating classroom results from A Taste of Genetics (**M6012**).

BLAST – Students will use BLAST to investigate the TAS2R38 protein sequences across various species and lineages to determine in which lineage did Type 2 taste receptor evolve, and use that information to construct a phylogenetic tree.

BSE – Don't Let the Cows Go Mad MiniLab

M6020

\$149



Regulations around cattle that exhibit Bovine Spongiform Encephalopathy (BSE), or Mad Cow disease, are strict and actions need to be swift when a potential contamination violation occurs. In this scenario a feed mill has self-reported a possible contamination. Students track down the source of contamination through species-specific PCR amplification and gel analysis, selecting which feed samples to test, and which controls to use.

Appropriate for high school biology students, especially honors and advanced placement, and college level biology.

Each MiniLab contains enough materials for 10 workstations, 2–3 students per workstation.

Materials include:

- | | |
|---|--|
| Ten 2% agarose GreenGel™ Cups | MiniOne® 5X Sample Loading Dye |
| FastTaq™ PCR MasterMix (2X) | MiniOne® DNA Marker |
| Primer set (forward and reverse) | One bag of 0.2 mL thin-wall PCR tubes |
| Three PCR control template DNA | One bag of 0.65 mL microcentrifuge tubes |
| DI Water | One bottle of 100 mL Tris-Borate-EDTA (TBE) buffer concentrate |
| Six DNA samples extracted from feed batches | Teacher's Guide |

Organize your PCR workstations

with the
MiniOne® T-Rack
Combo Pack (**M3143**)

(see pg. 11 for details)



MiniOne® Supplemental Reagent Packs

Designed for use in standardized biomedical curriculum programs that include a storyline, these supplemental reagent packs boost the activities as efficient hands-on labs. Each reagent pack includes enough materials for 10 workstations:

Sickle Cell Inheritance Reagent Pack


M3050

\$99

Students will run patient samples and compare against the controls to identify which family members are affected, carriers or normal for sickle cell anemia.

Materials include:

Ten 1.5% agarose GreenGel™ Cups	One bag of 0.65 mL microcentrifuge tubes
Eight DNA samples (three controls, five patients)	One bag of 2–200 µL micropipette tips
MiniOne® DNA Marker	Experimental protocol
One bottle of 100 mL TBE buffer concentrate	



Hypercholesterolemia Reagent Pack

M3051

\$99

Students will determine genotype for familial hypercholesterolemia for a family of five.

Materials include:

Ten 0.8% agarose GreenGel™ Cups	One bag of 0.65 mL microcentrifuge tubes
Seven DNA samples (two controls, five patients)	One bag of 2–200 µL micropipette tips
MiniOne® Universal DNA Marker	Experimental protocol
One bottle of 100 mL TBE buffer concentrate	



DNA Investigations Reagent Pack

M3052

\$120

Students will try to identify a missing person using restriction enzymes on DNA from two missing people and comparing to the restriction patterns from DNA from a skeleton.

Materials include:

Ten 1% agarose GreenGel™ Cups	One bag of 0.65 mL microcentrifuge tubes
Four DNA samples (Two Skeleton DNA pre-cut with EcoRV or HindIII, Undigested DNA from Missing Person's 1 and 2)	One bag of 2–200 µL micropipette tips
MiniOne® DNA Marker	MiniOne® 5X Sample Loading Dye
Two restriction enzymes	One bottle of 100 mL TBE buffer concentrate
One bag of 0.65 mL microcentrifuge tubes	Enzyme dilution buffer
	Experimental protocol



MiniOne® Gel Electrophoresis Starter Kit

M3200

\$75

This kit includes all supplies your class needs to get started with gel electrophoresis. Includes materials for making and running fifty MiniOne gels (1–2% agarose) - a great value! An exclusive MiniOne® DNA marker for analyzing a wide range of fragment sizes is also included.

Kit contents:

Item	Description
Agarose	10 grams, Electrophoresis Grade, Low EEO
TBE buffer concentrate	500 mL at 20X concentration
GelGreen™ DNA stain	50 µL at 10,000X stock
5X Sample Loading Dye	1 mL at 5X concentration, with Orange G and Xylene Cyanol tracking dyes
MiniOne® DNA Marker	500 µL for 50 loads, with 100, 300, 500, 1,000 and 2,000 bp bands in a ready-to-use format
1XTE buffer	2 mL, for DNA sample dilution
Reusable plastic gel cups	15 pieces, for making your own GreenGel™ Cups
Instruction manual	"How to Make MiniOne® Agarose Gels"



MiniOne® MiniLabs Alignment

Table of Courses/Grades and Recommended Kits

MiniLab Part #	MiniLab Name	Catalog Page	Middle School	Gen Bio	Honors Bio
ELECTROPHORESIS					
M3001	Electrophoresis 101	30	x	x	x
M3002	Gel Loading Practice Kit	28	x	x	x
M3003	PTC Genetics	31		x	x
M3004	DNA Fingerprinting	30	x	x	
M3005	CSI Forensics	31		x	
M3006	Foodborne Outbreak Investigation	34			x
M3007	Colorful Dye Electrophoresis	28	x	x	x
M3008	NGSS-Aligned Color Dyes and Gel Electrophoresis	35	x	x	
M3009	Candy Color Electrophoresis	29	x	x	
M3010	Hunting the Inheritance of Huntington's Disease	32			x
M3011	Determine the Genetics of a Ca\$H Cow	29	x	x	
M3012	PTC Inheritance and Graphical Analysis	32			
M6050	Restriction Digest Basics	33		x	x
M6053	Restriction Analysis of DNA	33			x
PCR and ELECTROPHORESIS					
M6001	PCR 101	36			x
M6010	A Taste of Genetics	38			x
M6005	PCR Cycle Number Analysis	37			
M6020	BSE - Don't Let the Cows Go Mad!	39			
REAGENT PACKS					
M3050	Sickle Cell Inheritance	40			
M3051	Hypercholesterolemia	40			
M3052	DNA Detectives Investigations	40			

AP Bio	Forensics	Health Science	Food Science	Agriculture	Environmental	Biotech	Biomed curriculum
x	x	x	x	x	x	x	
x	x	x	x	x	x	x	x
x		x				x	
	x				x		
	x						
x		x	x	x	x		
x						x	
			x				
x		x					
				x	x		
x						x	
x						x	
		x					x
		x					x
	x						x





MiniOne® GreenGel™ Cups

A ready to use, microwavable plastic cup with all ingredients to make one MiniOne gel.




- Each cup contains GelGreen™ DNA stain* mixed in a pre-weighed amount of agarose gel cubes
- One gel cup makes one MiniOne gel (10 mL in volume)
- Making agarose gels becomes easy, convenient, and fast
- 6 month shelf life


Procedure




Partially peel the plastic film



Microwave for 20 seconds



Pour into gel tray



Ready to run in 10–15 minutes

*Features of GelGreen™ DNA stain:

- Safety: a non-cytotoxic, non-mutagenic substitute for ethidium bromide (EtBr)
- Stability: stable at room temperature for long-term storage and microwavable
- Simplicity: binds to dsDNA, fluoresces when exposed to blue light enabling instant visualization of DNA bands

GreenGel™ Cups with GelGreen™ DNA stain for gel electrophoresis with DNA samples

Cat. No.	Buffer Type	Agarose Gel Concentration	Content	List Price
M3102TBE	TBE	1%	Ten GreenGel™ Cups with GelGreen™ DNA stain mixed in agarose gel cubes, and one bottle of 100 mL TBE buffer concentrate (Makes ten MiniOne® TBE gels)	\$24.00
M3142TBE	TBE	1.5%		\$24.00
M3103TBE	TBE	2%		\$24.00
M3123TBE	TBE	3%		\$25.00
M3141TBE	TBE	0.6%		\$24.00
M3140TBE	TBE	0.8%	\$24.00	
M3102TAE	TAE	1%	Ten GreenGel™ Cups with GelGreen™ DNA stain mixed in agarose gel cubes (Makes ten MiniOne® TAE gels)	\$21.00
M3142TAE	TAE	1.5%		\$21.00
M3103TAE	TAE	2%		\$21.00
M3123TAE	TAE	3%		\$22.00
M3141TAE	TAE	0.6%		\$21.00
M3140TAE	TAE	0.8%	\$21.00	

Agarose gel cups without DNA stain for gel electrophoresis with color dye samples

Cat. No.	Buffer Type	Agarose Gel Concentration	Content	List Price
M3151TBE	TBE	1%	Ten agarose gel cups with preweighed agarose gel cubes, and one bottle of 100 mL TBE buffer concentrate (Makes ten MiniOne® TBE gels)	\$20.00
M3151TAE	TAE	1%	Ten agarose gel cups with preweighed agarose gel cubes (Makes ten MiniOne® TAE gels)	\$18.00



DNA Size Markers


MiniOne® DNA Marker

M3104

\$35

A DNA size marker that consists of five double-stranded DNA fragments with sizes of 2K, 1K, 500, 300, and 100 base pairs (bp).

- 500 µL for 50 loads (10 µL per load)
- Suitable for 1% and 2% agarose gels
- All DNA bands will be well separated within 25 minutes




100 bp DNA Ladder

M3117

\$65

Consists of ten double-stranded DNA fragments, sizes ranging from 100 bp to 1,000 bp in exact 100 bp increments.

- 1,000 µL for 100 loads (10 µL per load)
- A size marker for most PCR products smaller than 1000 bp




MiniOne® Universal DNA Marker

M3144

\$35

A DNA size marker uniquely designed for fast band separation on agarose gels. It is composed of nine double-stranded DNA fragments with sizes of 10K, 6K, 3K, 2K, 1K, 800, 600, 400, and 200 base pairs (bp), reference band at 1 kb.

- 500 µL for 50 loads (10 µL per load)
- All DNA bands will be well separated within 25 minutes in a 1% agarose gel
- Suitable to be used as a size marker for most PCR products and recombinant plasmids and inserts




500 bp DNA Ladder

M3145

\$35

Consists of ten double-stranded DNA fragments, sizes ranging from 500 bp to 5,000 bp in exact 500 bp increments, reference band at 2.5 kb.

- 500 µL for 50 loads (10 µL per load)
- An ideal size marker for PCR products smaller than 5 kb




1 kb DNA Ladder

M3116

\$65

Consists of fifteen double-strand DNA fragments, size ranging from 1 kb to 15 kb in exact 1 kb increments, reference band at 5 kb.

- 1,000 µL for 100 loads (10 µL per load)
- Can be used as a size marker for restriction digestions of genomic DNA, large dsDNA fragment



All DNA markers are supplied in ready-to-load format with Xylene Cyanol FF and Orange G tracking dyes, stable for six months at room temperature.



Consumables and Plastics

Consumables and Plastics

Consumables and Plastics

TBE Buffer Concentrate (20X)
M3101TBE

\$18 (500 mL)

- Suitable for separation of smaller size DNA fragments or PCR products (≤ 2 kb)



TAE Buffer Concentrate (10X)
M3101TAE

\$18 (500 mL)

- Suitable for separation of larger size DNA fragments (5-20 kb), example: restriction digests of Lambda DNA



5X Sample Loading Dyes with Orange G and Xylene Cyanol
M3115

\$10 (10 mL)

5X Sample Loading Dyes with Orange G
M3119

\$10 (10 mL)

- Mix with DNA sample for easy loading
- Tracking dyes give a green color to samples which can be seen easily with the blue light on
- Other commonly used sample loading dyes are dark blue color which makes the DNA sample invisible when the blue light is on



FastTaq™ DNA Polymerase

A specially engineered Taq DNA polymerase with a very fast PCR extension rate at 100 bp/second. Also possesses moderate 3'-5' proofreading activity, making this enzyme well suited for high-throughput PCR.

MiniOne® FastTaq™ PCR MasterMix (2X)
M6201

\$90 (5 x 1 mL)

- The mastermix includes FastTaq™ DNA polymerase, dNTPs, Mg²⁺ ions and buffer
- Just add primers and template DNA to complete the reaction setup
- Sufficient for 500x10 uL reactions

MiniOne® FastTaq™ DNA Polymerase
M6202

\$90 (1,000 units at 5U/ μ L)

- FastTaq™ DNA Polymerase (5 U/ μ L): 0.2 mL
- 5X PCR Buffer, with Mg⁺⁺: 4 mL
- 25 mM MgSO₄: 1 mL



0.2 mL PCR Tubes
M6100

\$10

- Thin-walled with attached flat cap, optically clear, non-sterile, natural color
- Pack of 100 tubes
- Nuclease free



Taq DNA Polymerase

A regular Taq DNA polymerase suitable for a wide range of DNA assays with excellent yield and sensitivity. Routine PCR amplification of DNA templates up to 6 kb with a fast PCR extension rate at 1,000 bp/minute.

MiniOne® Taq PCR MasterMix (2X)
M6208

\$70 (5 mL)

- The mastermix includes Taq DNA polymerase, dNTPs, Mg²⁺ ions and buffer
- Just add primers and template DNA to complete the reaction setup



MiniOne® Taq DNA Polymerase
M6207

\$65 (1,000 units at 5U/ μ L)

- Enzyme is separately supplied with 10X PCR buffer and 25 mM MgSO₄

Micropipette Tips—Bulk package

M3112 1-10 μ L, pk of 250 tips
\$9.50

M3111 2-200 μ L, pk of 250 tips
\$9.50

M3134 2-200 μ L, pk of 1,000 tips
\$28.50

M3118 100-1,000 μ L, pk of 250 tips
\$9.50

- Fine tip with standardization marks
- Universal fit
- Autoclavable
- Non-sterile



PCR Reagents

You may find the following items handy when doing PCR:

MiniOne® dNTP Mix
M6203

\$55 (1 mL)

- A ready-to-use aqueous solution containing dATP, dCTP, dGTP and dTTP, each at a final concentration of 10 mM

Nuclease free water for PCR
M6204

\$8 (5 mL)



Nuclease free water
M6205

\$15 (50 mL)

- Molecular biology grade
- For DNA sample dilution or general use

1X Tris-EDTA (TE) Buffer, pH8.0
M6206

\$15 (50 mL)

- Molecular biology grade
- For DNA sample dilution

Micropipette Tips—Racked package

M3136 2-200 μ L, ten racks
\$65 \$55

- Tips packed in racks of 96 micropipette tips for easy use and storage
- Has all features of the bulk packaging version



GelGreen™ DNA Stain (10,000X concentration)
M3113

\$20 (50 μ L)

M3114

\$100 (500 μ L)

M3120

\$180 (1 mL)

M3121

\$850 (5 mL)

- Simple to use: just add 1 μ L per 10 mL agarose solution to make one MiniOne gel
- Safe: a non-cytotoxic, non-mutagenic, and environmentally safe substitute for ethidium bromide (EtBr)
- Stable at room temperature
- Microwavable



Agarose Electrophoresis grade, low EEO
M3105

\$15 (5 grams)

M3106

\$45 (25 grams)

M3106-100g

\$100 (100 grams)

M3106-500g

\$400 (500 grams)



Microcentrifuge Tubes
M3107

\$9.50 (0.65 mL, natural color)

M3109

\$9.50 (1.7 mL, natural color)

M3108

\$10 (0.65 mL, rainbow colors)

M3110

\$10 (1.7 mL, rainbow colors)

- Pack of 200 tubes
- Non-sterile



Visit our MiniOne® Resource Center

<https://theminione.com/minione-resource-center/>



Grant/Funding Resources

Grant templates and Donor's Choose hints



Instructional Videos and User Manuals

Equipment walk through, app set and GelCup prep



Free Classroom Activities

Free activities from MiniOne® Systems and our collaborators build on the skills and concepts learned in MiniOne® MiniLabs




MiniOne® Blog

What are teachers saying about MiniOne in their classroom?



General Information

Company information and education partners



MiniLab Teacher and Student Guides

Library of our MiniLab instructions



Professional Development

In-house teacher training



Conference Calendar

See if we're heading to a local conference in your area




Partnerships

Programs that provide outreach, loaners, teacher training and summer workshops



Index

42V power supply.....	7	Microcentrifuge tubes.....	48
Agar plates.....	28	Micropipette tips.....	49
Agarose.....	48	Micropipettes.....	12-15
Agarose gel cups.....	44, 45	MiniLabs.....	28-37
Amber filter.....	6, 17	MiniLabs alignment.....	42-43
A taste of genetics.....	43	MiniOne® Buffer Tank.....	7
Bioinformatics.....	38	MiniOne® Carriage Unit.....	6
BLAST.....	38	MiniOne® Centrifuges.....	16
Blue LED illuminator.....	17	MiniOne® Gel Casting System.....	8
BSE.....	39	MiniOne® Micropipettes.....	12
Buffers.....	48	MiniOne® PCR/Electrophoresis packages.....	24-27
Ca\$h cow.....	29	MiniOne® Photo Hood.....	6
Centrifuges.....	16	Mobile controller.....	23
Color dye labs.....	28, 30, 35	Multi-speed centrifuge.....	16
Combs.....	9	NCBI.....	38
Course alignment.....	42, 43	NGSS-aligned minilabs.....	35
Crime scene investigation.....	31	Nuclease free water.....	49
CSI.....	31	ONE Series™ Micropipettes.....	14-15
Deoxynucleotide triphosphates.....	49	PCR 101 minilab.....	36
DNA amplification.....	20, 22, 36	PCR 101 minilab and gel electrophoresis combo.....	36
DNA fingerprinting.....	30	PCR app.....	22
DNA investigations.....	40	PCR cycle number analysis.....	37
DNA ladders.....	46, 47	PCR mastermix.....	49
DNA markers.....	46	PCR polymerase.....	49
DNA stain.....	48	PCR reagents.....	49
dNTPs.....	49	PCR system.....	18-21
Electrophoresis 101.....	30	PCR tubes.....	49
Electrophoresis starter kit.....	41	PCR water.....	49
Electrophoresis system.....	4-6, 24-26	Photo hood.....	6, 17
Erlenmeyer flasks.....	11	Pipette stands.....	13, 15
FastTaq™ DNA polymerase.....	49	Pipettes.....	12-15
Food poisoning.....	34	Power supply.....	4, 7, 18
Foodborne outbreak investigation.....	34	PrepOne™ Sapphire.....	17
Forensics.....	31, 40	PTC genetics.....	31
Gel cups.....	44, 45	PTC inheritance.....	38
Gel electrophoresis starter kit.....	41	Ready-to-load DNA markers.....	46, 47
Gel loading practice.....	28	Restriction analysis of DNA.....	33
Gel tray platforms.....	7	Restriction digest basics.....	33
Gel trays.....	9	Sample loading dyes.....	48
GelGreen™ DNA stain.....	48	Sickle cell inheritance.....	40
Gel stain.....	48	Single speed centrifuge.....	16
Graphical analysis.....	32	Size markers.....	46, 47
GreenGel™ Cups.....	44, 45	Supplemental reagent packs.....	40
Hands-on labs.....	28-37	TAE buffer.....	48
Hardy-Weinberg.....	38	Taq DNA polymerase.....	49
Hypercholesterolemia.....	40	TBE buffer.....	48
Informatics.....	38	TE buffer.....	49
Lambda phage.....	36	The Winston™.....	2
Laminated practice pipette sheets.....	28	The Winston™ Platform.....	2
Laminated color dye cards.....	29	T-Rack™ Systems.....	10-11
Mad cow disease.....	39		



miniOne®

7738 Arjons Drive, San Diego, CA 92126 USA
1 858 684-3190 / 1 800 255-1777
1 858 684-3195 (fax)

 theminiOne.com

 info@theminiOne.com

©2020 Embi Tec. All rights reserved. Embi Tec and MiniOne are registered trademarks of C.C. IMEX.
FastTag, GreenGel, and PrepOne are trademarks of Embi Tec. GelGreen is a trademark of Biotium.
Bluetooth is a registered trademark of Bluetooth SIG, Inc. Contributing illustration from Science Lab Studios, Inc.



CAT-08-0320