


miniOne®

SYSTEMS

Product Catalog 2020 International





The Legacy of Winston Walker

Winston Walker dedicated his life to serving his country and finding solutions, whether it be on his aircraft or in industry after his extensive military service. His drive to make science education more accessible led to the development of the MiniOne® Electrophoresis and PCR Systems. His final contribution is The Winston viewer, which helps students of all ages answer various questions around molecules. His impact on science education is global, bringing equality to teachers and students worldwide.

MiniOne® The Winston™ M1050

See the GLOW! Use fluorescence to detect and see your samples glow in **colors** you won't believe!

Molecules are hard to see with the naked eye—they can be colorless, too small, or not abundant enough. Fluorescence to the rescue! Fluorescence helps you study the invisible. The Winston allows you and your students to investigate phenomena such as: "Does your sample have DNA? Prove it!", "Where is the DNA?", or "What color is chlorophyll"?

Place up to 4 samples on the base and cover with the MiniOne® Photo Hood to reveal the GLOW! Each unit includes:

- One MiniOne® The Winston™ Platform
- One built-in rechargeable battery
- One MiniOne® Photo Hood
- USB charging cable

As a part of the community that Winston helped bring together, we invite you to submit what you discover with The Winston and we will post those protocols online to share with others. For these community-submitted activities check out our website:

<https://www.theminione.com/the-winston>

To submit an activity please email us at: info@minione.com

The Winston will make you **GLOW!**

MiniOne® The Winston™ Platform M1051

Includes the base component and charging cable. Great add-on if you already have a MiniOne Photo Hood.



Table of Contents

Electrophoresis Systems.....	4
T-Rack System.....	10
Micropipettes.....	12
Microcentrifuges.....	16
PrepOne™ Sapphire and Photo Hood	17
PCR System.....	18
PCR and Electrophoresis Packages	24
MiniLabs, Electrophoresis.....	28
MiniLabs, PCR.....	36
MiniLab Course Alignment.....	42
GreenGel™ Cups	44
DNA Size Markers	46
Consumables and Plastics.....	48
Resources	50
Index	51

Schedule your order to ship when you need it!

Simply let us know future ship dates when you place your order.

Ordering Information

Orders must be placed by:

Email europa@theminione.com

Telephone +39 02 35980827

MiniOne® Electrophoresis System

M1000-EU, M1000-UK


Integrated, real-time electrophoresis system for running DNA labs in one classroom period.



Includes one each of the following items:

(See pages 6 to 9 for complete product details)

- MiniOne® Carriage with blue LED illumination and magnet-activated safety switch
- Gel tank with graphite electrodes
- 42V power supply, 100–240V input
- Amber photo hood for real-time viewing and capturing gel image
- Casting system with casting stand, two gel trays, two reversible combs for six and nine wells, and lid
- **FREE** 2–20 μ L variable volume micropipette (\$59 value!)
- Validation kit: Two GreenGel™ Cups, three DNA samples and TBE buffer concentrate (one per order)



Example results from
Electrophoresis System validation kit
~20 min. run time

MiniOne® Electrophoresis Classroom Package of 10 Systems

M1010-EU, M1010-UK

One classroom package is for a class of ten student groups, two to three students per group.

Includes:


- Ten sets of MiniOne® Electrophoresis System (**M1000-EU or M1000-UK**)
- Ten **FREE** 2–20 μ L variable volume micropipette
- Validation kit: Two GreenGel™ Cups, three DNA samples and TBE buffer concentrate (one per order)



MiniOne® Carriage M2007

Housing and control unit with sealed LED lights and power controls.


- + and – electrical contacts for carbon electrodes of gel tank
- Design ensures correct orientation fit
- Two rows of blue LEDs illuminate the gel from the sides
- Choose from two intensities—low light for loading, bright light for DNA viewing
- Built-in, magnet-activated on/off power switch



MiniOne® Photo Hood M2005

Molded photo hood for 360° viewing with ventilation slots and black imaging platform.

- Four embedded magnets activate power carriage
- Amber filter blocks blue light and transmits green light showing bright fluorescent DNA bands against a dark background
- Ventilation slots prevent condensation inside the hood



NEW!

MiniOne® Photo Hood Phone Platform M2017 (Set of 5)

Molded adapter allows you to place a phone on top of the MiniOne Photo Hood to easily take videos, time lapses, or still images of your electrophoresis run. (Photo Hood not included)

- Latches into the photo hood vents
- Can support mobile devices up to 6.2" (L) x 3.0" (W)" long
- Hands-free documentation



MiniOne® 42V Power Supply M2006-EU, M2006-UK

Input voltage: 100–240 V, 50/60 Hz, 0.5A


- Output voltage: 42V, 0.19A
- European or UK plug
- Certifications: CE, FCC, cULus, RoHS compliant



MiniOne® Gel Tank M2001

Molded polycarbonate tank with graphite electrodes.

- Three ridges on walls of tank to guide correct orientation of gel tray
- Durable 1 cm graphite electrodes replace standard fragile platinum wires for easy cleaning
- Black and silver gray gel platform provides contrasting background for imaging gels




MiniOne® Gel Tray Platforms

M2014 Black gel tray platform (10 per pack)

M2015 Silver gray gel tray platform (10 per pack)

Slots on the sides of the platform ensure that it only fits into the tank in the correct orientation.




Use with DNA separation. Black color plastic plate printed with a fluorescent green ruler on the edge, texture on well area.

Use with color dyes separation. Silver gray color plastic plate with texture on well area.

MiniOne® Gel Casting System M2002

A compact and self-contained unit with two clear gel trays and two gel combs that fit inside a casting stand with lid.

- Molded casting stand eliminates leakage—no tape needed
- Gel tray guides create a one-way fit for the gel trays when casting, ensuring the proper positioning of the wells for use in the gel tank
- Dual reversible comb with one beveled side, 1.5 mm max. thickness, six and nine wells, and white lid for overnight storage of pre-poured gels
- Beveled comb for easy comb removal and sample loading
- White cover blocks light so you can keep the gels you pour ahead of time safe from light exposure



MiniOne® Gel Trays

M2013 (10 per pack)

M2013-Bulk (50 per pack)

New gel tray design, clear acrylic trays with one orientation fit feature.

- Use in MiniOne tank and casting stand



NEW!

MiniOne® Casting Stand Cover

M2018 (5 per pack)

Fits MiniOne Gel Casting Stand.

- Blocks light for advanced gel prep




MiniOne® Gel Combs

M2004 (Set of 2)

Dual reversible gel comb with one beveled side.

- 1.5 mm max. thickness, 6+6 / 9+9 wells



NEW!

Colorful MiniOne® T-Rack™ System

The modular MiniOne T-Racks are designed to give you the most flexibility with how you set up your lab stations.

MiniOne® T-Rack™ Microtube Rack - 1.5 and 2.0 mL

M3181 (Set of 5 racks)

Five assorted colors, holds 3 x 4 of 1.5 or 2.0 mL microtubes per rack. (Tubes not included)



Colors may vary




Capacity and configuration guide

MiniOne® T-Rack™ PCR Tube Rack - 0.2 and 0.5 mL

M3180 (Set of 5 racks)

Five assorted colors, each rack holds 2 x 5 of 0.2 mL PCR tubes and 3 x 4 of 0.5 mL microtubes. (Tubes not included)



Colors may vary



Capacity and configuration guide

MiniOne® T-Rack™ Micropipette Tip Rack and Cover

M3139 (Set of 5 racks)

Five assorted rack colors with clear lids, pre-racked with 2–200 µL Universal tips, low binding, 36 tips/rack.



Colors may vary



Capacity and configuration guide

MiniOne® T-Rack™ System Tray

M3182 (Set of 5 trays)

Trays for holding up to three T-Rack modules per platform for easy workstation organization.



MiniOne® T-Rack Combo Pack

M3143 (Set of 5)

Assorted colors (5), complete set of racks for 0.2 mL/0.5 mL PCR tubes, 1.5/2.0 mL microtubes, and 2–200 µL universal tips, with convenient tray to hold up to 3 racks at a time. Combo pack contains one set each of M3139, M3180, M3181, M3182, to make 5 complete sets.

**NEW!**

MiniOne® Mini Erlenmeyer Flasks

M2019 (Set of 5)

Safer than traditional glass, this affordable set of five polypropylene 170 mL Erlenmeyer flasks is ideal for the biology classroom. They're light, rugged, autoclavable, and feature a pour spout, making them perfect for storing and pouring buffer. Minimize spills and breakage in the lab with MiniOne Mini Erlenmeyer Flasks. Non-microwavable, do not freeze.

- 170 mL (6 oz) graduated every 25 mL (1 oz)
- Polypropylene
- Autoclavable (*remove sticker before autoclaving*)



MiniOne® Micropipettes

M2008, M2010, M2011, M2012
Best Classroom Pipette at the Best Price!

One micropipette per pack; comes with inspection certificate, calibration tool, operation manual and sample pipette tips.



Features:

- Adjustable volume micropipette with durable and reliable quality
- Ergonomic design provides comfortable operation for small or large hands
- Two 'stops' on the plunger to allow for accurate and complete sample dispensation
- Ultra-affordable for educational use
- Can be self-calibrated

MiniOne® Micropipette Set

M2016



This set of three lab-quality micropipettes comes protected in a case constructed of durable, easy-to-clean, waterproof-coated Lycra, ensuring your pipette set stays safe, organized, and clean.

Includes one of each variable volume micropipettes:

M2008 2–20 µL

M2010 20–200 µL

M2011 100–1000 µL

MiniOne® Micropipette Stand

M2021

Keep your bench organized!
 Holds up to nine MiniOne Micropipettes.



Cat. No.	MiniOne® Micropipette Volume Range	Accuracy	Repeatability
M2008	2–20 µL (H20)	± 2.5-1.0%	≤ 1.50–0.30%
M2010	20–200 µL (H200)	± 1.8-0.6%	≤ 0.50–0.15%
M2011	100–1000 µL (H1000)	± 1.5-0.6%	≤ 0.30–0.15%
M2012	1–10 µL (H10)	± 2.5-1.0%	≤ 1.50–0.40%
M2016	Set of 3 MiniOne® Micropipettes, 2–20 µL, 20–200 µL and 100–1000 µL adjustable volume micropipettes (one each) in a protective carrying case		
M2021	Pipette stand for the MiniOne® micropipettes; hold 9 micropipettes		

All micropipettes are manufactured according to the standards ISO13485. QC inspection and test are complying with standard ISO8655.

Micropipette Tips - See page 49 for specifications

M3112	1–10 µL, pk of 250 tips 1–10 µL, pk of 250 tips
M3111	2–200 µL, pk of 250 tips
M3134	2–200 µL, pk of 1,000 tips
M3118	100–1,000 µL, pk of 250 tips
M3139	Micropipette Tip Rack and Cover, set of 5 assorted colors <ul style="list-style-type: none"> • 36 x 2–200 µL Universal tip per rack • See page 10 for details

ONE Series™ Micropipettes

EA-1001 to EA-1006

One micropipette per pack; comes with inspection certificate, calibration tool, operation manual and sample pipette tips.



Features:

- Adjustable volume micropipette combines all the features of accuracy, reliability, and ergonomic design in one
- Great value for scientific research and educational uses
- Autoclavable
- Can be self-calibrated

ONE Series™ Micropipette Stand

EA-1011

Keep your bench organized!
Holds up to nine ONE Series Micropipettes.



Cat. No.	ONE Series™ Micropipette Volume Range	Accuracy	Repeatability
EA-1001	0.1–2 µL	± 12.0–15%	≤ 6.00–0.70%
EA-1002	2–20 µL	± 2.5–1.0%	≤ 1.50–0.30%
EA-1003	20–200 µL	± 1.8–0.6%	≤ 0.50–0.15%
EA-1004	10–100 µL	± 1.8–0.8%	≤ 0.50–0.15%
EA-1005	100–1000 µL	± 1.5–0.6%	≤ 0.30–0.15%
EA-1006	0.5–10 µL	± 2.5–1.0%	≤ 1.50–0.40%
EA-1010	Set of 3 ONE Series Micropipettes: 2–20 µL, 20–200 µL and 100–1000 µL adjustable volume micropipettes and MiniOne pipet holder		
EA-1011	Pipette Stand for the ONE Series Micropipettes; holds 9 micropipettes		
EA-1014	Set of 4 ONE Series Micropipettes: 0.1–2 µL, 2–20 µL, 20–200 µL and 100–1000 µL adjustable volume micropipettes (one each), w/Holder		

All micropipettes are manufactured according to the standards ISO13485. QC inspection and test are complying with standard ISO8655.

Micropipette Tips - See page 49 for specifications	
M3112	1–10 µL, pk of 250 tips 1–10 µL, pk of 250 tips
M3111	2–200 µL, pk of 250 tips
M3134	2–200 µL, pk of 1,000 tips
M3118	100–1,000 µL, pk of 250 tips
M3136	2–200 µL, ten racks <ul style="list-style-type: none"> • Tips packed in racks of 96 micropipette tips for easy use and storage • Has all features of the bulk packaging version

MiniOne® Centrifuges

Portable mini centrifuges ideal for bringing small droplets to the bottom of tubes, for micro-filtrations, or basic separations. Brushless motor for low noise level. CE marked.

MiniOne® Multi Speed Centrifuge

M2031

Features:

- Adjustable speed from 1,000 to 10,000 RPM
- One combi-rotor with positions for varied sizes of microcentrifuge tubes and PCR tubes—no need to change rotors



Specifications

Input Voltage	100-240V AC, 50-60Hz
Rotor Capacity	2 x 8 x 0.2 mL PCR tubes and 6 x 1.5/2.0 mL and 6 x 0.5 mL microcentrifuge tubes
Speed Range	Adjustable speed: 1,000 to 10,000 RPM, max 3,200 x g
Timer	15 sec to 99 minutes or continuous operation
Dimensions	20 x 16 x 13 cm
Weight	1.05 kg

MiniOne® Single Speed Centrifuge

M2032

Features:

- Easy operation: one button to control ON or OFF
- Fixed speed at 10K RPM providing 4,800 x g (RCF) when using the microcentrifuge tube rotor



Specifications

Input Voltage	100-240V AC, 50-60Hz
Rotor Capacity	PCR tube rotor: 2 x 8 x 0.2 mL PCR tubes, Microcentrifuge tube rotor: 6 x 1.5/2.0 mL or 6 x 0.5 mL with adaptors or 6 x 0.2 mL with adaptors
Speed Range	Fixed speed, 10,000 RPM
Dimensions	17 x 15 x 13 cm
Weight	0.95 kg

PrepOne™ Sapphire and Photo Hood M5000

A blue LED illuminator for smart devices to take gel images with ease. Adapt your current electrophoresis system for safe blue light illumination and non-toxic GelGreen™ stain.

Package includes:

- One PrepOne™ Sapphire Blue LED Illuminator
- One photo hood for smart devices (assembly required)
- One amber filter

Photo hood


- Compact and collapsible design turns your bench into an instant dark room.
- Wide top platform holds most smart devices
- Height: 15.4 cm, provides a suitable focal length for most phone cameras
- Base DIM: 13.8 (W) x 13.8 (L) cm

Amber filter

- Cuts out excess blue light to enhance the contrast of the DNA signal
- DIM: 15.2 (W) x 15.2 (L) x 0.3 (H) cm

PrepOne™ Sapphire

- Side illuminating blue light allows direct visualization of DNA bands
- 4x4 blue LEDs last over 50,000 hours
- Two built-in light intensities
- DIM: 21.4 (W) x 16.3 (L) x 1.5 (H) cm





MiniOne® PCR System

M4000-EU, M4000-UK


Teach and Do PCR Labs in 45 minutes!

Fast PCR thermal cycler controlled by App on Bluetooth® LE mobile device (see page 22). Amplify DNA in one classroom period.

- One MiniOne thermal cycler for PCR with 16 sample wells
- One 12V, 100W power supply, 100–240 VAC
- One **FREE** validation kit
- Temperature range **4°C (below ambient) to 99°C**
- Ability to pause run to evaluate cycle number



Validation kit provides enough reagents for three PCR runs, each amplifying three different size fragments, then visualizing the results with MiniOne® Electrophoresis



Example results from PCR System validation kit ~ 20 min. amplification and ~ 20 min. run time

Lane 1: MiniOne DNA marker (see page 34)
Lane 2: PCR fragment 1
Lane 3: PCR fragment 2
Lane 4: PCR fragment 3
Lane 5: negative control

Features:

- Peltier cooling technology and custom algorithm drives **fast thermal** cycling
- **Mobile App** for programming and monitoring via **Bluetooth® LE**
- Indicator lights on front show status of the run
- Fully compatible with standard reagents, consumables, and protocols

Benefits:

- Complete a PCR protocol in a single class period
- Intuitive programming interface for students
- Keep your samples at 4°C when done



Specifications

Sample capacity	16 x 0.2 mL standard PCR tubes
Temperature range	4°–99°C
Heated lid	Yes, with safety interlock
Communications	Bluetooth® Low Energy wireless technology
Software	Graphical programming interface; Real-time protocol monitoring
Weight	1.9 lb (860 g) approx.
Dimensions	12 x 12 x 12 cm (4.7 x 4.7 x 4.7 in) approx.
Operating voltage	100–240 VAC



MiniOne® PCR System (continued)

16 sample capacity

Perfect for two to four student groups or a small class.

For larger classes, add more systems to maintain student exposure and participation.



Active heating and cooling

The Peltier element drives rapid temperature transitions and faster protocols. Keeps your samples at 4°C at the end of the run.

Specialized algorithm

Custom control algorithm is the key for fast cycling and precise thermal control.

Visual indicators

Bright LEDs show the state of the machine and progress of the protocol.



Compact

The MiniOne® PCR System takes up minimal space on your bench, stores easily, and can be moved wherever it's needed.



Safe

Fully enclosed system for safe operation in the classroom.




Heated Lid

No condensation or messy mineral oil. Safety switch keeps protocol from starting until the lid is closed.

How fast is it?

A typical fast protocol - MiniOne® PCR System is done in **50%** of the time.

Step	Temperature	Duration	Cycles
Denaturation	94°C	5 sec	20 cycles
Annealing	54°C	5 sec	
Extension	72°C	5 sec	



MiniOne® PCR App

Intuitive, student-centered app for programming and monitoring your PCR protocols.

Features:

- Interactive screens lead students through each step of setting up the PCR protocol
- **Pause** feature stops the protocol at the end of the extension step for convenient cycle number analysis
- Graphical output **displays real-time temperature data** that can be saved or emailed at the end of the run
- **Linked protocol** to automatically run consecutive protocols

1 Bluetooth® LE connection



View and connect to available PCR Systems with a wireless Bluetooth® LE connection. No cables or wifi needed!

2 Setup



Create a new protocol, access your last protocol or browse the library of saved protocols from the Setup menu.

3 Protocol




Select constant temperature mode for restriction digests, transformation and DNA extractions, or PCR mode for DNA amplification.

4 Constant Temperature



Set up constant temperature protocols to incubate samples at 0-99°C using a simple, intuitive interface.

5 PCR



Graphic PCR programming interface allows students to enter time and temperatures directly on the PCR overview graph.

6 Cycling



The run screen shows a continuously updated, real-time readout of the temperature, current cycle, current step, and elapsed time.



iPhone
(iOS 9 or
higher)



Android phone
(Lollipop or
higher)



iPad/Mini iPad



MiniOne® PCR



Android controller



MacBook or PC Laptop*



Chromebook
(Bluetooth 4.1 or above)

*Requires MiniOne® Bluetooth® LE Dongle **M4060**



Android Mobile Controller M4050

Android mobile controller with MiniOne® PCR App pre-loaded.

- 7-inch, full color graphic display for a modern, student-centered interface
- Bluetooth® LE connectivity for programming and monitoring your MiniOne® PCR System



MiniOne® PCR with Android Mobile Controller

M4001

- One MiniOne® PCR System with one validation kit - **CE** marked
- One Android Mobile Controller with MiniOne® PCR App installed



MiniOne® PCR/Electrophoresis Package I

M4011

- One MiniOne® PCR System with one validation kit - **CE** marked
- One MiniOne® Electrophoresis System - **CE** marked
 - One **FREE** 2–20 μ L variable volume micropipette

For a group of two to three students



Package II

M4012

- One MiniOne® PCR System with one validation kit - **CE** marked
- Two MiniOne® Electrophoresis Systems - **CE** marked
 - Two **FREE** 2–20 μ L variable volume micropipettes

For two groups of students



Package III

M4026

- Two MiniOne® PCR Systems with two validation kits - **CE** marked
- Six MiniOne® Electrophoresis Systems - **CE** marked
 - Six **FREE** 2–20 μ L variable volume micropipettes

For six groups of students



MiniOne® PCR/Electrophoresis Package IV

M4039

Three MiniOne® PCR Systems - **CE** marked

- Two PCR validation kits
- Nine MiniOne® Electrophoresis Systems - **CE** marked
- Nine **FREE** 2–20 µL variable volume micropipettes

For nine groups
of students



MiniOne® Bioscience Classroom Starter Package V

M4040

This bioscience classroom starter package includes the essential equipment for teaching hands-on molecular biology concepts using DNA amplification and separation in the classroom. The starter package is ideal for nine groups of two to three students, and includes:

- Three MiniOne® PCR Systems - **CE** marked
- Two PCR validation kits
- Nine MiniOne® Electrophoresis Systems - **CE** marked
- Three MiniOne® Single Speed Microcentrifuges - **CE** marked
- One MiniOne® Micropipette Set with Case, includes one of each: 2–20, 20–200, 100–1000 µL variable volume micropipettes
- Nine **FREE** 2–20 µL variable volume micropipettes

For nine groups
of students

Includes
micropipette set
for teacher prep!



MiniOne® MiniLabs

Our hands-on MiniLabs are a fun and engaging series of modules that take students from mastery of basic biotech skills, through popular applications of electrophoresis in forensics, DNA fingerprinting, and human genetics, and finally, to a challenging, real-world investigation of a foodborne outbreak.

Gel Loading Practice MiniLab

M3002

Teach your students how to read, adjust, and use a micropipette. Practice pipetting and loading samples into the wells of a real gel before handling valuable DNA samples.

For new users.

Each MiniLab contains enough materials for 10 workstations, 2–3 students per workstation.

Materials include:

Twenty precast agar plates with 3 rows of 8 wells per row	One bag of 1.7 mL microcentrifuge tubes
Twenty adhesive practice pipetting sheets	One bag of 2–200 μ L micropipette tips
Two tubes of dye samples (10 mL each)	Teacher's Guide



Laminated Practice Pipette Cards

M2022 (Set of 20)

Reusable practice pipette cards from our Gel Loading Practice MiniLab. Lamination makes it easy to clean so you are ready for your back-to-back classes.



Colorful Dye Electrophoresis MiniLab

M3007


This fun and colorful lab introduces students to gel electrophoresis principles, including the basics of electricity and macromolecules. Predict how molecular size and electrical charge affect a molecule's migration in a separation matrix.

Ideal for ages 11–14.

Each MiniLab contains enough materials for 10 workstations, 2–3 students per workstation.

Materials include:

Ten 1% agarose gel cups	One bag of 0.65 mL microcentrifuge tubes
Nine color dye samples	One bag of 2–200 μ L micropipette tips
One bottle of 100 mL Tris-Borate-EDTA (TBE) buffer concentrate	Teacher's Guide



Each MiniLab is an all-inclusive science lab kit. The ready-to-pour gel cups greatly simplify steps to make agarose gels. The TBE buffer concentrate and easy-to-load DNA/color dye samples assure good results, minimize teacher prep work, and maximize student success in the lab.

Candy Color Electrophoresis MiniLab

M3009

Investigate the phenomenon of food dyes using candy to illustrate the effect of mass and charge during electrophoresis.

Ideal for ages 11–14.

Each MiniLab contains enough materials for 10 workstations, 2–3 students per workstation.

Materials include:

Ten 1% agarose gel cups	Dye extraction buffer
One bottle of 100 mL Tris-Borate-EDTA (TBE) buffer concentrate	Candies in six colors
One bag of 1.7 mL microcentrifuge tubes	Ten 10-well dye extraction trays
One bag of 2–200 μ L micropipette tips	Teacher's Guide



Laminated Gel Annotation and Photo Template

M2023 (Set of 10)

Place your gel on these laminated cards make the color dyes POP after your electrophoresis run! Dual sided for 6- or 9-well dye electrophoresis activities.



Determining the Genetics of a Ca\$H Cow MiniLab

M3011

Using an understanding behind cheese production, genetic inheritance, and gel electrophoresis, students will determine genotype of two bulls and three cows, and recommend which combination a dairy farmer should purchase to produce more valuable offspring.

Ideal for ages 11–17.

Each MiniLab contains enough materials for 10 workstations, 2–3 students per workstation.

Materials include:

Eight color dye samples	One bag of 0.65 mL microcentrifuge tubes
Ten 1.5 % agarose gel cups	One bag of 2–200 μ L micropipette tips
One bottle of 100 mL Tris-Borate-EDTA (TBE) buffer concentrate	Teacher's Guide



Electrophoresis 101 MiniLab

M3001


Students are introduced to the principles of gel electrophoresis by separating colorful dyes and DNA samples on an agarose gel. Challenge your students' analytical and mathematical skills as they construct a standard curve to determine the sizes of unknown DNA fragments.

Ideal for ages 14–17.

Each MiniLab contains enough materials for 10 workstations, 2–3 students per workstation.

Materials include:

- Ten 1% agarose GreenGel™ Cups
- Three color dye samples and three DNA samples
- One bottle of 100 mL Tris-Borate-EDTA (TBE) buffer concentrate
- One bag of 0.65 mL microcentrifuge tubes
- One bag of 2–200 µL micropipette tips
- Teacher's Guide



DNA Fingerprinting MiniLab

M3004

How is DNA used to trace the history and heritage of an individual? Students help scientists identify the father of a baby humpback whale using DNA fingerprinting technology. Engage your students with a real-world application of genetics as they analyze a complex array of DNA bands to arrive at a logical solution.

Appropriate for ages 14–17.

Each MiniLab contains enough materials for 10 workstations, 2–3 students per workstation.

Materials include:

- Ten 1% agarose GreenGel™ Cups
- Five DNA samples
- One bottle of 100 mL Tris-Borate-EDTA (TBE) buffer concentrate
- One bag of 0.65 mL microcentrifuge tubes
- One bag of 2–200 µL micropipette tips
- Teacher's Guide



CSI Forensics MiniLab

M3005


Explore a crime scene investigation in the classroom. Learn an exciting real-world application of gel electrophoresis and the statistical principles of human genetic identification. Students will logically integrate multiple lines of evidence, including fingerprints, hair samples, and DNA fingerprinting to connect an individual to a crime scene and solve the mystery of "Who Killed Dr. Ward?"

Appropriate for ages 14–17.

Each MiniLab contains enough materials for 10 workstations, 2–3 students per workstation.

Materials include:

- Ten 1% agarose GreenGel™ Cups
- Five DNA samples
- One bottle of 100 mL Tris-Borate-EDTA (TBE) buffer concentrate
- One bag of 0.65 mL microcentrifuge tubes
- One bag of 2–200 µL micropipette tips
- Teacher's Guide



PTC Genetics MiniLab


M3003

Solve a genetic mystery using gel electrophoresis. Students develop a hypothesis about the inheritance of a trait in a family, then test their hypotheses by running restriction fragments on a gel, analyzing a Punnett square, and constructing a family tree.

Appropriate for ages 14–17.


Each MiniLab contains enough materials for 10 workstations, 2–3 students per workstation. Materials include:

- Ten 2% agarose GreenGel™ Cups
- Six pre-digested DNA samples
- Forty pieces of PTC tasting and taste control papers
- One bottle of 100 mL TBE buffer concentrate
- One bag of 0.65 mL microcentrifuge tubes
- One bag of 2–200 µL micropipette tips
- Teacher's Guide



PTC Inheritance and Graphical Analysis MiniLab

M3012



Explore Mendelian genetic inheritance, use Punnett Squares to make a claim, and see your evidence with DNA electrophoresis. This lab also features DNA fragment size analysis and demystifies why the standard curve is plotted as Log-Y!

Appropriate for ages 14-17 and advanced courses.

Each MiniLab contains enough materials for 10 workstations, 2-3 students per workstation.

Materials include:

Ten 2% agarose GreenGel™ Cups
Six pre-digested DNA samples
MiniOne® DNA Marker
Forty pieces of PTC tasting and taste control papers

One bottle of 100 mL Tris-Borate-EDTA (TBE) buffer concentrate
One bag of 0.65 mL microcentrifuge tubes
One bag of 2-200 µL micropipette tips
Teacher's Guide

Hunting the Inheritance of Huntington's Disease MiniLab

M3010



In this lab, students will examine family history to construct a pedigree and will assess molecular data to make predictions about inheritance of the disease in fraternal twins. Students will then perform gel electrophoresis to confirm genotype.

Appropriate for ages 14-17 and advanced courses.

Each MiniLab contains enough materials for 10 workstations, 2-3 students per workstation.

Materials include:

Ten 2% agarose GreenGel™ Cups
Four Ready-to-Load DNA samples
100 bp DNA Ladder
One bottle of 100 mL Tris-Borate-EDTA (TBE) buffer concentrate
One bag of 0.65 mL microcentrifuge tubes
One bag of 2-200 µL micropipette tips
Teacher's Guide

Restriction Digest Basics MiniLab

M6050



Cells have mechanisms for cutting long strands of nucleic acid into shorter strands—a type of molecular scissors. There are several reasons why cells need to cut their DNA or RNA. In this lab students will explore what restriction enzymes do and determine electrophoresis fragment sizes by comparing bands to a molecular weight standard.

Appropriate for ages 14-17.

Each MiniLab contains enough materials for 10 workstations, 2-3 students per workstation.

Materials include:

Ten 1.5% agarose GreenGel™ Cups
Four DNA samples (three pre-digested, one undigested)
MiniOne® Universal DNA Marker

One bottle of 100 mL Tris-Borate-EDTA (TBE) buffer concentrate
One bag of 0.65 mL microcentrifuge tubes
One bag of 2-200 µL micropipette tips
Teacher's Guide

Restriction Analysis of DNA MiniLab

M6053



"Molecular scissors" for cutting DNA or RNA are used by cells for many reasons, including defense against an invading host, or genetic recombination. These restriction enzymes allow researchers to study smaller pieces of DNA more thoroughly. In this restriction digestion lab, students will explore what restriction enzymes do, perform a single and double digest of a synthetic piece of DNA, predict fragment sizes and compare to fragment size of the actual digested DNA run on an agarose gel.

Appropriate for ages 14-17 and advanced courses.

Each MiniLab contains enough materials for 10 workstations, 2-3 students per workstation.

Materials include:

Ten 1.5% agarose GreenGel™ Cups
Undigested DNA Sample
Four pre-digested DNA controls
MiniOne® Universal DNA marker
Enzyme dilution buffer
Two restriction enzymes
MiniOne® 5X Sample Loading Dye

One bottle of 100 mL Tris-Borate-EDTA (TBE) buffer concentrate
One bag of 0.65 mL microcentrifuge tubes
One bag of 2-200 µL micropipette tips
One bag of 0.2 mL PCR tubes
Teacher's Guide

Foodborne Outbreak Investigation MiniLab M3006

This MiniLab is a student-driven discovery process based on a real *Shigella* outbreak in 2000. Students use scientific reasoning and forensic science principles to analyze epidemiological data, develop a hypothesis, and test their hypothesis with gel electrophoresis. They will tabulate data expressed as text to systematically analyze the case and evaluate experimental approaches used by their team and others in the class. Along the way they will develop an understanding of foodborne outbreaks and foodborne illness, topics that are frequently in the news and relevant to students' lives.

Appropriate for ages 14-17 and advanced courses.

Each MiniLab contains enough materials for 10 workstations, 2-3 students per workstation.

Materials include:

- Ten 1% agarose GreenGel™ Cups
- Eleven DNA samples
- One bottle of 100 mL Tris-Borate-EDTA (TBE) buffer concentrate
- Two bags of 0.65 mL microcentrifuge tubes
- One bag of 2-200 µL micropipette tips
- Teacher's Guide
- A two-part detailed guide with background information, step-by-step procedures and worksheets for students
- Quiz question bank and answers



NGSS-Aligned Color Dyes and Gel Electrophoresis MiniLab M3008



A comprehensive, 5E inquiry, week-long lesson plan that introduces students to separation science in a variety of contexts. Student-driven inquiry activities impart the scientific background needed to understand gel electrophoresis. The curriculum culminates in a gel electrophoresis experiment using colorful dyes. The curriculum is provided as a downloadable PDF manual which includes comprehensive background readings, student worksheets, and teacher support materials. Lab materials and a PowerPoint presentation are also included.

Ideal for ages 11-14.

Each MiniLab contains enough materials for 10 workstations, 2-3 students per workstation.

Materials include:

- Ten 1% agarose gel cups
- Six color dye samples
- One set of filter papers and green food coloring for exploring paper chromatography
- One bottle of 100 mL Tris-Borate-EDTA (TBE) buffer concentrate
- One bag of 0.65 mL microcentrifuge tubes
- One bag of 2-200 µL micropipette tips
- One copy of the curriculum



Start

10 minutes
run time
→



Finish

Bundle and Save!

M3008 3-Pack Reagents for 30 workstations, 1 copy of the curriculum	\$337
M3008 5-Pack – Reagents for 50 workstations, 1 copy of the curriculum	\$475

	1	3	7	9	10	11	2	4	5	6	8	12
Potato Chips		X	X	X		X	X	X				
French Onion Dip	X	X	X	X	X	X	X	X		X	X	X
Tortilla Chips		X	X	X								
Salsa	X			X	X			X			X	X
Guacamole				X			X	X	X	X		
Homemade Potato Salad	X	X	X	X	X	X	X	X				
Burgers						X	X	X				
Hot Dogs	X	X	X	X	X	X	X	X				
Five Layer Bean Dip			X									
Buttermilk	X	X	X	X	X	X	X	X				

PCR 101 MiniLab: Amplification from the Lambda Phage Genome

M6001



In this hands-on PCR MiniLab, students use polymerase chain reaction (PCR) to amplify three segments of the Lambda phage genome. They will look at sequence data, predict the fragment sizes of the PCR products, then compare their predictions to the PCR products they amplify and run on an agarose gel. Complete amplification in 17 minutes with the MiniOne® PCR system and the FastTaq™ Master Mix.

Appropriate for ages 14-17 and advanced courses.

Each MiniLab contains enough materials for 10 workstations, 2-3 students per workstation.

Materials include:

Ten 2% agarose GreenGel™ Cups
FastTaq™ PCR MasterMix (2X)
Three primer sets, both forward and reverse primers included in each set
Lambda phage genomic DNA
Sterile nuclease-free water
MiniOne® DNA marker

MiniOne® 5X Sample Loading Dye
One bag of 0.2 mL thin-wall PCR tubes
One bag of 0.65 mL microcentrifuge tubes
One bottle of 100 mL Tris-Borate-EDTA (TBE) buffer concentrate
Teacher's Guide

Bundle and Save!


	M6001	M6002	M6004	M6003
M6002 - PCR 101 MiniLab, PCR Reagents (no electrophoresis reagents)	✓	✓	✓	✓
M3103TBE - Ten 2% agarose GreenGel™ Cups with TBE buffer concentrate (see page 45)	✓			✓
M3136 - 10 racks of micropipette tips (2-200 µL) (see page 49)			✓	✓

COMPLETE PACKAGE!

PCR Cycle Number Analysis

MiniLab

M6005




Visualize the power of exponential growth with PCR! Students will set up PCR reactions and analyze the products after a variable number of cycles. Students will estimate the minimum number of cycles needed to detect a PCR product on an agarose gel, and visualize and appreciate exponential growth.

Appropriate for ages 14-17 and advanced courses.

Each MiniLab contains enough materials for 10 workstations, 2-3 students per workstation.

Materials include:

Ten 1% agarose GreenGel™ Cups
MiniOne® DNA marker
MiniOne® 5X Sample Loading Dye
FastTaq™ PCR MasterMix (2X)
One primer set, forward and reverse
Lambda phage genomic DNA
One bag of 0.65 mL microcentrifuge tubes
One bag of 0.2 mL thin-wall PCR tubes
One bottle of 100 mL TBE buffer concentrate
Teacher's Guide



Need some racked tips? or PCR tube racks?



Try the

MiniOne® T-Rack Micropipette Tip Rack and Cover (**M3139**) or the MiniOne T-Rack™ PCR Tube Rack - 0.2 and 0.5 mL (**M3180**) (see pg. 10 for details)

A Taste of Genetics MiniLab: Extract and Amplify the PTC Gene M6010



This hands-on MiniLab introduces students to the science of human genetic variation through DNA extraction, PCR amplification, restriction digest, and analysis of the TAS2R38 taster gene of their own DNA, and compare genotype to phenotype.

Appropriate for ages 14-17 and advanced courses.

Each MiniLab contains enough materials for 10 workstations, 2–3 students per workstation.

Materials include:

DNA extraction solution	One bag of 0.65 mL microcentrifuge tubes
Forward and reverse primers for PTC genes	Forty pieces of PTC taste paper and taste control papers
Taq polymerase master mix (2X)	Ten 2% agarose GreenGel™ Cups
HaeIII restriction enzyme	One bottle of 100 mL Tris-Borate-EDTA (TBE) buffer concentrate
Restriction enzyme dilution buffer	Two grams table salt
MiniOne® 5X Sample Loading Dye	Teacher's Guide
MiniOne® DNA Marker	
One bag of 0.2mL thin-walled PCR tubes	

Bundle and Save!

	M6010	M6012	M6013
M6012 - A Taste of Genetics MiniLab, DNA Extraction and PCR Reagents (no electrophoresis reagents)	✓	✓	✓
M3103TBE - Ten 2% agarose GreenGel™ Cups with TBE buffer concentrate (see page 45)	✓		✓
M3136 - 10 racks of micropipette tips (2–200 µL) (see page 49)			✓

COMPLETE PACKAGE!

Extension Activities

available at <https://theminione.com/free-classroom-activities/>

Introduction to NCBI Bioinformatics – Students will learn to navigate NCBI, explore the different types of information available, and apply it to look deeper at the TAS2R38 gene, including sequence, chromosomal location, and associated phenotypes and variants.

Hardy-Weinberg – Students will create a simulation of allele frequencies in a population using principles from the Hardy-Weinberg model and a TI-84 Plus graphing calculator or computer including evaluating classroom results from A Taste of Genetics (**M6012**).

BLAST – Students will use BLAST to investigate the TAS2R38 protein sequences across various species and lineages to determine in which lineage did Type 2 taste receptor evolve, and use that information to construct a phylogenetic tree.

BSE– Don't Let the Cows Go Mad MiniLab M6020



Regulations around cattle that exhibit Bovine Spongiform Encephalopathy (BSE), or Mad Cow disease, are strict and actions need to be swift when a potential contamination violation occurs. In this scenario a feed mill has self-reported a possible contamination. Students track down the source of contamination through species-specific PCR amplification and gel analysis, selecting which feed samples to test, and which controls to use.

Appropriate for ages 14-17 and advanced courses.

Each MiniLab contains enough materials for 10 workstations, 2–3 students per workstation.

Materials include:

Ten 2% agarose GreenGel™ Cups	MiniOne® 5X Sample Loading Dye
FastTaq™ PCR MasterMix (2X)	MiniOne® DNA Marker
Primer set (forward and reverse)	One bag of 0.2 mL thin-wall PCR tubes
Three PCR control template DNA	One bag of 0.65 mL microcentrifuge tubes
DI Water	One bottle of 100 mL Tris-Borate-EDTA (TBE) buffer concentrate
Six DNA samples extracted from feed batches	Teacher's Guide

Organize your PCR workstations

with the
MiniOne® T-Rack
Combo Pack (**M3143**)

(see pg. 11 for details)



MiniOne® Supplemental Reagent Packs

Designed for use in standardized biomedical curriculum programs that include a storyline, these supplemental reagent packs boost the activities as efficient hands-on labs. Each reagent pack includes enough materials for 10 workstations:

Sickle Cell Inheritance Reagent Pack

M3050



Students will run patient samples and compare against the controls to identify which family members are affected, carriers or normal for sickle cell anemia.

Materials include:

Ten 1.5% agarose GreenGel™ Cups	One bag of 0.65 mL microcentrifuge tubes
Eight DNA samples (three controls, five patients)	One bag of 2–200 µL micropipette tips
MiniOne® DNA Marker	Experimental protocol
One bottle of 100 mL TBE buffer concentrate	

Hypercholesterolemia Reagent Pack

M3051



Students will determine genotype for familial hypercholesterolemia for a family of five.

Materials include:

Ten 0.8% agarose GreenGel™ Cups	One bag of 0.65 mL microcentrifuge tubes
Seven DNA samples (two controls, five patients)	One bag of 2–200 µL micropipette tips
MiniOne® Universal DNA Marker	Experimental protocol
One bottle of 100 mL TBE buffer concentrate	

DNA Investigations Reagent Pack

M3052



Students will try to identify a missing person using restriction enzymes on DNA from two missing people and comparing to the restriction patterns from DNA from a skeleton.

Materials include:

Ten 1% agarose GreenGel™ Cups	One bag of 0.65 mL microcentrifuge tubes
Four DNA samples (Two Skeleton DNA pre-cut with EcoRV or HindIII, Undigested DNA from Missing Person's 1 and 2)	One bag of 2–200 µL micropipette tips
MiniOne® DNA Marker	MiniOne® 5X Sample Loading Dye
Two restriction enzymes	One bottle of 100 mL TBE buffer concentrate
One bag of 0.65 mL microcentrifuge tubes	Enzyme dilution buffer
	Experimental protocol

MiniOne® Gel Electrophoresis Starter Kit

M3200

This kit includes all supplies your class needs to get started with gel electrophoresis. Includes materials for making and running fifty MiniOne gels (1–2% agarose) - a great value! An exclusive MiniOne® DNA marker for analyzing a wide range of fragment sizes is also included.

Kit contents:

Item	Description
Agarose	10 grams, Electrophoresis Grade, Low EEO
TBE buffer concentrate	500 mL at 20X concentration
GelGreen™ DNA stain	50 µL at 10,000X stock
5X Sample Loading Dye	1 mL at 5X concentration, with Orange G and Xylene Cyanol tracking dyes
MiniOne® DNA Marker	500 µL for 50 loads, with 100, 300, 500, 1,000 and 2,000 bp bands in a ready-to-use format
1XTE buffer	2 mL, for DNA sample dilution
Reusable plastic gel cups	15 pieces, for making your own GreenGel™ Cups
Instruction manual	"How to Make MiniOne® Agarose Gels"





MiniOne® MiniLabs Alignment

Table of Courses/Grades and Recommended Kits

MiniLab Part #	MiniLab Name	Catalog Page	Life Sciences
ELECTROPHORESIS		Ages 11-14	
M3001	Electrophoresis 101	30	x
M3002	Gel Loading Practice Kit	28	x
M3003	PTC Genetics	31	
M3004	DNA Fingerprinting	30	x
M3005	CSI Forensics	31	
M3006	Foodborne Outbreak Investigation	34	
M3007	Colorful Dye Electrophoresis	28	x
M3008	NGSS-Aligned Color Dyes and Gel Electrophoresis	35	x
M3009	Candy Color Electrophoresis	29	x
M3010	Hunting the Inheritance of Huntington's Disease	32	
M3011	Determine the Genetics of a Ca\$H Cow	29	x
M3012	PTC Inheritance and Graphical Analysis	32	
M6050	Restriction Digest Basics	33	
M6053	Restriction Analysis of DNA	33	
PCR and ELECTROPHORESIS			
M6001	PCR 101	36	
M6010	A Taste of Genetics	38	
M6005	PCR Cycle Number Analysis	37	
M6020	BSE - Don't Let the Cows Go Mad!	39	
REAGENT PACKS			
M3050	Sickle Cell Inheritance	40	
M3051	Hypercholesterolemia	40	
M3052	DNA Detectives Investigations	40	

Gen Bio	Honors Bio	AP/IB Biology	Forensics	Health Science	Food Science	Agriculture	Environmental	Biotech
Ages 14-17+								
x	x	x	x	x	x	x	x	x
x	x	x	x	x	x	x	x	x
x	x	x		x				x
x			x				x	
x			x					
	x	x		x	x	x	x	
x	x	x						x
x								
x					x			
	x	x		x				
x						x	x	
		x						x
x	x	x						
	x	x						x
	x	x						x
	x	x		x				x
		x				x		x
		x			x	x		x
				x				
				x				
			x					



MiniOne® GreenGel™ Cups

A ready to use, microwavable plastic cup with all ingredients to make one MiniOne gel.




- Each cup contains GelGreen™ DNA stain* mixed in a pre-weighted amount of agarose gel cubes
- One gel cup makes one MiniOne gel (10 mL in volume)
- Making agarose gels becomes easy, convenient, and fast
- 6 month shelf life


Procedure




Partially peel the plastic film



Microwave for 20 seconds



Pour into gel tray



Ready to run in 10–15 minutes

*Features of GelGreen™ DNA stain:

- Safety: a non-cytotoxic, non-mutagenic substitute for ethidium bromide (EtBr)
- Stability: stable at room temperature for long-term storage and microwavable
- Simplicity: binds to dsDNA, fluoresces when exposed to blue light enabling instant visualization of DNA bands

GreenGel™ Cups with GelGreen™ DNA stain for gel electrophoresis with DNA samples

Cat. No.	Buffer Type	Agarose Gel Concentration	Content
M3102TBE	TBE	1%	Ten GreenGel™ Cups with GelGreen™ DNA stain mixed in agarose gel cubes, and one bottle of 100 mL TBE buffer concentrate (Makes ten MiniOne® TBE gels)
M3142TBE	TBE	1.5%	
M3103TBE	TBE	2%	
M3123TBE	TBE	3%	
M3141TBE	TBE	0.6%	
M3140TBE	TBE	0.8%	Ten GreenGel™ Cups with GelGreen™ DNA stain mixed in agarose gel cubes (Makes ten MiniOne® TAE gels)
M3102TAE	TAE	1%	
M3142TAE	TAE	1.5%	
M3103TAE	TAE	2%	
M3123TAE	TAE	3%	
M3141TAE	TAE	0.6%	
M3140TAE	TAE	0.8%	

Agarose gel cups without DNA stain for gel electrophoresis with color dye samples


Cat. No.	Buffer Type	Agarose Gel Concentration	Content
M3151TBE	TBE	1%	Ten agarose gel cups with preweighed agarose gel cubes, and one bottle of 100 mL TBE buffer concentrate (Makes ten MiniOne® TBE gels)
M3151TAE	TAE	1%	Ten agarose gel cups with preweighed agarose gel cubes (Makes ten MiniOne® TAE gels)

DNA Size Markers

MiniOne® DNA Marker M3104

A DNA size marker that consists of five double-stranded DNA fragments with sizes of 2K, 1K, 500, 300, and 100 base pairs (bp).


- 500 µL for 50 loads (10 µL per load)
- Suitable for 1% and 2% agarose gels
- All DNA bands will be well separated within 25 minutes



MiniOne® Universal DNA Marker M3144

A DNA size marker uniquely designed for fast band separation on agarose gels. It is composed of nine double-stranded DNA fragments 3K, 2K, 1K, 800, 600, 400, and 200 base pairs (bp), reference band at 1 kb.


- 500 µL for 50 loads (10 µL per load)
- All DNA bands will be well separated within 25 minutes in a 1% agarose gel
- Suitable to be used as a size marker for most PCR products and recombinant plasmids and inserts



1 kb DNA Ladder M3116

Consists of fifteen double-strand DNA fragments, size ranging from 1 kb to 15 kb in exact 1 kb increments, reference band at 5 kb.


- 1,000 µL for 100 loads (10 µL per load)
- Can be used as a size marker for restriction digestions of genomic DNA, large dsDNA fragments



100 bp DNA Ladder M3117

Consists of ten double-stranded DNA fragments, sizes ranging from 100 bp to 1,000 bp in exact 100 bp increments.


- 1,000 µL for 100 loads (10 µL per load)
- A size marker for most PCR products smaller than 1000 bp



500 bp DNA Ladder M3145

Consists of ten double-stranded DNA fragments, sizes ranging from 500 bp to 5,000 bp in exact 500 bp increments, reference band at 2.5 kb.

- 500 µL for 50 loads (10 µL per load)
- An ideal size marker for PCR products smaller than 5 kb



All DNA markers are supplied in ready-to-load format with Xylene Cyanol FF and Orange G tracking dyes, stable for six months at room temperature.



Consumables and Plastics

TBE Buffer Concentrate (20X)

M3101TBE
(500 mL)

- Suitable for separation of smaller size DNA fragments or PCR products (≤ 2 kb)



TAE Buffer Concentrate (10X)

M3101TAE
(500 mL)

- Suitable for separation of larger size DNA fragments (5–20 kb), example: restriction digests of Lambda DNA



5X Sample Loading Dyes with Orange G and Xylene Cyanol

M3115
(10 mL)

5X Sample Loading Dyes with Orange G

M3119
(10 mL)

- Mix with DNA sample for easy loading
- Tracking dyes give a green color to samples which can be seen easily with the blue light on
- Other commonly used sample loading dyes are dark blue color which makes the DNA sample invisible when the blue light is on



GelGreen™ DNA Stain (10,000X concentration)

M3113
(50 μ L)

M3114
(500 μ L)

M3120
(1 mL)

M3121
(5 mL)

- Simple to use: just add 1 μ L per 10 mL agarose solution to make one MiniOne gel
- Safe: a non-cytotoxic, non-mutagenic, and environmentally safe substitute for ethidium bromide (EtBr)
- Stable at room temperature
- Microwavable



Agarose Electrophoresis grade, low EEO

M3105
(5 grams)

M3106
(25 grams)

M3106-100g
(100 grams)

M3106-500g
(500 grams)



Microcentrifuge Tubes

M3107
(0.65 mL, natural color)

M3109
(1.7 mL, natural color)

M3108
(0.65 mL, rainbow colors)

M3110
(1.7 mL, rainbow colors)

- Pack of 200 tubes
- Non-sterile



FastTaq™ DNA Polymerase

A specially engineered Taq DNA polymerase with a very fast PCR extension rate at 100 bp/second. Also possesses moderate 3'–5' proofreading activity, making this enzyme well suited for high-throughput PCR.

MiniOne® FastTaq™ PCR MasterMix (2X)

M6201
(5 x 1 mL)

- The mastermix includes FastTaq™ DNA polymerase, dNTPs, Mg^{2+} ions and buffer
- Just add primers and template DNA to complete the reaction setup
- Sufficient for 500x10 μ L reactions

MiniOne® FastTaq™ DNA Polymerase

M6202
(1,000 units at 5U/ μ L)

- FastTaq™ DNA Polymerase (5 U/ μ L): 0.2 mL
- 5X PCR Buffer, with Mg^{++} : 4 mL
- 25 mM $MgSO_4$: 1 mL



0.2 mL PCR Tubes

M6100

- Thin-walled with attached flat cap, optically clear, non-sterile, natural color
- Pack of 100 tubes
- Nuclease free



Taq DNA Polymerase

A regular Taq DNA polymerase suitable for a wide range of DNA assays with excellent yield and sensitivity. Routine PCR amplification of DNA templates up to 6 kb with a fast PCR extension rate at 1,000 bp/minute.

MiniOne® Taq PCR MasterMix (2X)

M6208
(5 mL)

- The mastermix includes Taq DNA polymerase, dNTPs, Mg^{2+} ions and buffer
- Just add primers and template DNA to complete the reaction setup



MiniOne® Taq DNA Polymerase

M6207
(1,000 units at 5U/ μ L)

- Enzyme is separately supplied with 10X PCR buffer and 25 mM $MgSO_4$

Micropipette Tips–Bulk package

M3112
1–10 μ L, pk of 250 tips

M3111
2–200 μ L, pk of 250 tips

M3134
2–200 μ L, pk of 1,000 tips

M3118
100–1,000 μ L, pk of 250 tips

- Fine tip with standardization marks
- Universal fit
- Autoclavable
- Non-sterile



PCR Reagents

You may find the following items handy when doing PCR:

MiniOne® dNTP Mix

M6203

(1 mL)

- A ready-to-use aqueous solution containing dATP, dCTP, dGTP and dTTP, each at a final concentration of 10 mM

Nuclease free water for PCR

M6204

(5 mL)

Nuclease free water

M6205

(50 mL)

- Molecular biology grade
- For DNA sample dilution or general use



1X Tris-EDTA (TE) Buffer, pH8.0

M6206

(50 mL)

- Molecular biology grade
- For DNA sample dilution

Micropipette Tips–Racked package

M3136 2–200 μ L, ten racks

- Tips packed in racks of 96 micropipette tips for easy use and storage
- Has all features of the bulk packaging version





Visit our MiniOne® Resource Center

<https://theminione.com/minione-resource-center/>



Grant/Funding Resources

Grant templates and Donor's Choose hints



Instructional Videos and User Manuals

Equipment walk through, app set and GelCup prep



Free Classroom Activities

Free activities from MiniOne® Systems and our collaborators build on the skills and concepts learned in MiniOne® MiniLabs




MiniOne® Blog

What are teachers saying about MiniOne in their classroom?



General Information

Company information and education partners



MiniLab Teacher and Student Guides

Library of our MiniLab instructions



Professional Development

In-house teacher training



Conference Calendar

See if we're heading to a local conference in your area




Partnerships

Programs that provide outreach, loaners, teacher training and summer workshops



Index

42V power supply.....	7	Microcentrifuge tubes.....	48
Agar plates.....	28	Micropipette tips.....	49
Agarose.....	48	Micropipettes.....	12-15
Agarose gel cups.....	44, 45	MiniLabs.....	28-37
Amber filter.....	6, 17	MiniLabs alignment.....	42-43
AP Biology.....	43	MiniOne® Buffer Tank.....	7
A taste of genetics.....	38	MiniOne® Carriage Unit.....	6
Bioinformatics.....	38	MiniOne® Centrifuges.....	16
BLAST.....	38	MiniOne® Gel Casting System.....	8
Blue LED illuminator.....	17	MiniOne® Micropipettes.....	12
BSE.....	39	MiniOne® PCR/Electrophoresis packages.....	24-27
Buffers.....	48	MiniOne® Photo Hood.....	6
Ca\$h cow.....	29	Mobile controller.....	23
Centrifuges.....	16	Multi-speed centrifuge.....	16
Color dye labs.....	28, 30, 35	NCBI.....	38
Combs.....	9	NGSS-aligned minilabs.....	35
Course alignment.....	42, 43	Nuclease free water.....	49
Crime scene investigation.....	31	ONE Series™ Micropipettes.....	14-15
CSI.....	31	PCR 101 minilab.....	36
Deoxynucleotide triphosphates.....	49	PCR 101 minilab and gel electrophoresis combo.....	36
DNA amplification.....	20, 22, 36	PCR app.....	22
DNA fingerprinting.....	30	PCR cycle number analysis.....	37
DNA investigations.....	40	PCR mastermix.....	49
DNA ladders.....	46, 47	PCR polymerase.....	49
DNA markers.....	46	PCR reagents.....	49
DNA stain.....	48	PCR system.....	18-21
dNTPs.....	49	PCR tubes.....	49
Electrophoresis 101.....	30	PCR water.....	49
Electrophoresis starter kit.....	41	Photo hood.....	6, 17
Electrophoresis system.....	4-6, 24-26	Pipette stands.....	13, 15
Erlenmeyer flasks.....	11	Pipettes.....	12-15
FastTaq™ DNA polymerase.....	49	Power supply.....	4, 7, 18
Food poisoning.....	34	PrepOne™ Sapphire.....	17
Foodborne outbreak investigation.....	34	PTC genetics.....	31
Forensics.....	31, 40	PTC inheritance.....	38
Gel cups.....	44, 45	Ready-to-load DNA markers.....	46, 47
Gel electrophoresis starter kit.....	41	Restriction analysis of DNA.....	33
Gel loading practice.....	28	Restriction digest basics.....	33
Gel tray platforms.....	7	Sample loading dyes.....	48
Gel trays.....	9	Sickle cell inheritance.....	40
GelGreen™ DNA stain.....	48	Single speed centrifuge.....	16
Gel stain.....	48	Size markers.....	46, 47
Graphical analysis.....	32	Supplemental reagent packs.....	40
GreenGel™ Cups.....	44, 45	TAE buffer.....	48
Hands-on labs.....	28-37	Taq DNA polymerase.....	49
Hardy-Weinberg.....	38	TBE buffer.....	48
Hypercholesterolemia.....	40	TE buffer.....	49
Informatics.....	38	The Winston™.....	2
Lambda phage.....	36	The Winston™ Platform.....	2
Laminated practice pipette sheets.....	28	T-Rack™ Systems.....	10-11
Laminated color dye cards.....	29		
Mad cow disease.....	39		



miniOne®

+39 02 35980827

 theminione.com

 europe@theminione.com

©2020 Embi Tec. All rights reserved. Embi Tec and MiniOne are registered trademarks of C.C. IMEX.
FastTaq, GreenGel, and PrepOne are trademarks of Embi Tec. GelGreen is a trademark of Biotium.
Bluetooth is a registered trademark of Bluetooth SIG, Inc. Contributing illustration from Science Lab Studios, Inc.



CAT-09-0520 -INTL